



PRODUCT GUIDE 2019



Because Time Is Everything

WHY ZIRC

2

COLOR METHOD™

4 - 19

Why the Color Method™?	4
Color Chart	4
Zobe	5
Tub Kits	7
Tub Accessories	8
Trays Covers	10
Multi-Mod Racks	12
Cassettes	13
Steri-Containers Mats	15
Bur & Endo Guards/Blocks	16
E-Z ID Rings/Tape	19

ISOLATION

20 - 27

Isolation360™ Power Pack	20
Airway Armor®	22
Pink Petal®	23
mr.thirsty® one-step	24
Insti-Dam® SE Cushion Lingua-Fix®	26

MOUTH MIRRORS

28 - 33

Crystal HD® Mirrors	28
mirror magic® Anti-Fog Solution	32

OPERATORY/MISC

34 - 44

CAD/CAM Block Locker™	34
-----------------------	----

MAKING DENTAL TEAMS HAPPIER AND MORE EFFICIENT



Zirc History - 1967

James Campion, our CEO, founded Zirc in 1967 and his daughter Nicolle, who joined the team in 2000, is the company's president. In 2017 we celebrated 50 years of offering our customers the highest quality products at the best value. Dental technology and processes have seen many changes over the last 50 years, but at Zirc, our commitment to providing time-saving, innovative dental products continues to be our top priority. We look forward to serving our customers for years to come!



Color Method™

Zirc has a brilliant solution for organizing the dental practice by the method of color-coding. We have rebranded our color-coding system to the Color Method™, to better educate the concept of color-coding in the practice. The method gives step-by-step instructions to achieve an efficient, productive, and compliant work place...all while remaining simple and stress-free.



Made in the USA

Our headquarters, located in Buffalo, Minnesota, houses our engineering, tooling, manufacturing, distributing, and customer service operations. We're proud to say that our state-of-the-art, high-quality products are made in the U.S.A.



Antimicrobial

Antimicrobial compounds are intended to protect parts against staining, odors, and degradation due to fungi, mold, yeasts, mildew, algae, and other microbes. Zirc incorporates antimicrobial compounds during the manufacturing process. Unlike a coating, this prevents fading or wearing away, giving you long-lasting protection. The incorporation of antimicrobial compounds in our plastic products is not intended to protect users or others against food-borne or disease-causing bacteria. For more information, please reference EPA Registration Notice 2000-1.




Zirc Perks

Members receive FREE samples, money-saving promotions, tips to get and stay organized, product evaluation invites, and more! Sign up to become a Zirc Perks Member via www.zirc.com or call us at 800-682-6636 (available in U.S.A. only).



What's New?



ASPIRATION & INGESTION PREVENTION

Airway Armor®

Airway Armor® prevents foreign objects and debris from going down the throat during procedures, giving yourself a peace of mind and your patients the protection they deserve. It is easy to place, tear resistant, and comfortable to the patient, all while giving you full mouth accessibility.
See page 22.

Continuous, hands-free suction!



GIVE A HAND TO HYGIENE

Pink Petal®

Provides continuous, comfortable, hands-free suction with no additional attachments required. By placing buccal side, the Pink Petal® confidently holds the saliva ejector right where you need it.
See page 23.

Find the products that fit your needs!



STRUGGLING WITH ISOLATION?

Isolation360™ Power Pack

Zirc's Isolation360™ has you covered! The pack contains 7 different isolation devices, allowing you to become familiar with each one. Also included inside the pack, is an Ideal usage chart with product recommendations for different procedures.

Once you use all the products, order more of what works best for your practice!
See page 20.

CONTROL BLOCK CHAOS

CAD/CAM Block Locker™

Our CAD/CAM Block Locker™ is a stylish, modern cabinet designed to help organize the chaos of milling blocks. Utilize the online template to label each internal organizer tray by shade, color and opacity...creating a seamless way to check inventory, stock and keep blocks organized with ease.
See page 34.



For convenience, safety and efficiency!



Safe-Lok® Cover Tall

Maintain compliance with OSHA and CDC with our locking tray cover! Now offered in a tall version, you will be able to fit our Performance Tool Organizer® and almost any stainless steel cassette!
See page 11.

COMING SOON!

Pouch-it™

Zirc will be releasing a line of custom-fit pouches that will pair with our instrument holders, eliminating the need of guessing what size to get!
See page 13.



**Organization is the goal.
Color is the method.
Zirc is your solution.**

ORGANIZE WITH EASE
Why the Color Method™?

COLOR METHOD™
THE KEY TO ORGANIZATION



Zirc has the tools you need to manage your materials and easily store, access, and clean your instruments. Staff, and patients alike, will appreciate a streamlined, stress-free day, from set-up to sterilization. Our exclusive Color Method™ will help you get organized, increase efficiency and eliminate leaving chairside to search for materials.

Color coding has been proven to save time for each procedure, allowing you to see more patients in a day.

Simplified Practice

Once you implement the Color Method™, your workflow will become a systematic routine between the sterilization and treatment room, creating maximum efficiency!

Benefits

- Set-ups are organized, eliminating the need to find missing materials or instruments during a procedure
- Makes re-stocking dental materials more efficient by decreasing the amount of inventory, and reducing the amount of expired materials
- Simplified training of new staff
- Reduces the risk of an exposure incident by containing dirty instruments in covered a container, satisfying the CDC and OSHA standards
- Visually pleasing to patients versus stainless steel

Color Chart

With 16 product colors, there is guaranteed to be a combination that matches your practice!



Want help implementing the Color Method™?

For more information, visit us online to download our Color Method™ brochure.



COLOR METHOD™



web: www.zirc.com | p: +1 763.682.6636 | f: +1 763.682.6604

AM Antimicrobial

136 Autoclavable

Washer Disinfect

S Surface Disinfect



Syringe Composite Tub



Capsule Composite Tub



Endo Tub



Crown Prep Tub



Hygiene Tub



Implant Tub

SMALL CHANGES Big Results

Cluttered drawers = stress and chaos. Stay efficient with color-coded organization.

Start by removing everything from your operatory drawers. Place in a large box and set aside. Restock your operatory drawers with staple supplies such as: bibs, saliva ejectors, cotton products, etc. (basically, anything you use in EVERY procedure).

Here's where the Tubs come into play. You need a minimum of one Tub for each procedure. Each Tub should be a different color so that you can quickly identify procedures by color. Stock each Tub with all materials you will need for that procedure.

Store ALL Tubs in a central area (such as sterilization). Take the Tub into the treatment room when a procedure is performed. Pull out what you need at chairside. Disinfect and re-store.

Reduce your stress and have fun with color organization.





MAKE IT SIMPLE WITH AN All-In-One Kit

Our kits are designed to store and organize materials and keep your inventory under control. In addition, they are easily portable for transportation between procedures and your centralized location for sterilization, disinfecting, and restocking.


Complete Tub

Includes: Procedure Tub, Tub Liner, Tub Dividers, Cup/Cover Single, Cup/Cover Double, Divided Slide Tray, Safe-Lok Cover for Tub®.

Dimensions: 32cm x 30cm x 10.6cm

(20Z455) Available in all colors

AM S



Use: Any Procedure

Composite Kits

Store and organize composite materials with our convenient Composite Kits. When stored in a Tub, contents become portable and can be transported to a central sterilization area for restocking. Organizational accessories help you keep better inventory control with Clip-On Holders that allow you to separate composite for use during procedures. The Long Tub Cup and Slide Tray provide additional organizational space. The Shade Label Sheet provides self-stick labels preprinted with composite shade.

Syringe Composite Kit

Syringe holder allows you to store up to 40 syringe tubes. Kit includes: Syringe 20-Unit Stand, 1 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Syringe Clip-On Holder, Shade Label Sheet 2pk.

(20Z489)

AM S



Capsule Composite Kit

Long Tub Cups provide a total of 20 compartments for capsule storage. Includes: Tub Liner, 4 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Capsule Clip-On Holder, Shade Label Sheet 2pk.

(20Z479)

AM S



Build Your Own Tub



Procedure Tub

Allows you to “build” your own setup.
Dimensions: 31.12cm x 27.62cm x 6.98cm
(20Z463) Available in all colors



Safe-Lok® Tub Cover

Locking cover for securing materials.
Keep materials protected during storage and transportation.
Dimensions: 32.39cm x 30cm x 3.5cm
(20Z459)



Tub Liner

Foam liner that sits in bottom of Tub.
Tub not included. Replace as needed.
Dimensions: 24.45cm x 22.54cm x .32cm
(20Z462)



Tub Box

Designed to hold/organize “hard-to-find” supplies.
Non-slip bottom for secure placement.
Dimensions: 21.27cm x 5.87cm x 2.86cm
(20Z474) Gray Only



Tub Insert

10 compartments for easy organization.
Sits on rim of Tub, allowing for storage underneath.
Dimensions: 27.3cm x 24.6cm x 2.54cm
(20Z465) Gray Only



Tub Dividers

Creates up to 24 compartments.
May be cut or trimmed to create larger compartments.
Dimensions: 32.39cm x 30cm x 3.5cm
(20Z461) Gray Only



Use color-coded tape on Tub Cups and Tub Box to identify procedural contents. (Learn more on page 19)



Tub Cups with Covers

Fill Tub Cups with brushes, clamps, cotton pellets, prophylaxis angle tips, prophylaxis paste, pumice, syringe tips, implant tools, etc. Can be used in procedure tubs, drawers, on trays, etc. for small, easy-to-lose items. Covers are transparent.



Single Tub Cup with Cover

Dimensions: 5.27cm x 4.13cm x 2.86cm
(20Z471) Gray Only

Double Tub Cup with Cover

Dimensions: 8.6cm x 5.72cm x 2.86cm
(20Z472) Gray Only

Long Tub Cup with Cover

Dimensions: 21.27cm x 5.72cm x 2.86cm
(20Z473) Gray Only



Divided Slide Tray
 Dimensions: 24.77cm x 13.65cm x 2.7cm
 (20Z466) Clear

Flat Slide Tray
 Dimensions: 24.77cm x 13.65cm x 2.7cm
 (20Z467) Clear




Shade Label Sheets 2pk
 Preprinted labels for marking Tub materials by procedure. Shade Labels can be used with the Syringe Composite and Capsule Composite Tub Organizers.
 (20Z463-LBL)

Slide Trays

Sits on rim of Tub and allows for storage underneath.



Check out our Implant and Restorative Organizers on page 44.



Syringe 20-Unit Stand

Stands up for easy access to composite.
 Lays flat to fit inside Tub or drawer.
 Holds up to 20 syringe tubes.
 Has space for 2 Syringe 10-Unit Add-Ons (Stand can hold up to 40 syringes).
 Dimensions: 22.22cm x 12.06cm x 3.96cm
 (20Z487)



Syringe 10-Unit Add-on

An Add-on for the Syringe 20-Unit Stand.
 Gives additional space for syringe storage.
 Holds 10 syringe tubes.
 Dimensions: 12.74cm x 5.87cm x 1.9cm
 (20Z488)



Syringe 20-Unit Stand

The Syringe 10-Unit Add-on may be added for storage of up to 40 syringe tubes.



Syringe Clip-On Holder

Allows quick access to composite during procedure.
 Conveniently clips on any side of Tub, B-Lok Tray or Capsule Composite.
 Holds 15 syringe tubes and 3 applicator tips.
 Dimensions: 20.32cm x 6.19cm x 1.9cm
 (20Z486)



Capsule Clip-On Holder

Allows quick access to composite during procedure.
 Conveniently clips on any side of Tub or B-Lok Tray.
 Holds 8 capsules, 4 syringe tubes, and capsule/gun dispenser.
 Dimensions: 20.32cm x 5.08cm x 1.42cm
 (20Z476)



Trays & Covers



Are you in compliance?

Trays

Our exclusive Safe-Lok® tray cover provides an excellent method of infection control and safety that allows for two benefits other trays do not:

- 1** Prevents any contaminated items from falling off of the used tray during transportation from operator to central sterilization. The contaminated items stay enclosed until you are ready for the cleaning and sterilization process.
- 2** Set up your procedure trays ahead of time with all instrumentation needed. The Safe-Lok® cover securely locks to the tray, not only keeping items clean, but also allowing for safe transportation to chairside. The cover is removed once the patient is in the chair and the procedure begins.

OSHA

Immediately or as soon as possible after use, contaminated reusable sharps shall be placed in puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. (To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.)

CDC

Dirty Instruments must be “contained” prior to transporting them. Containment can be defined as **instrument cassettes** or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.

The Safe-Lok® Tray Cover helps prevent cross contamination and provide safe transportation back to the sterilization area.





B-Lok Flat Tray

Wide edges for secure grip.
Non-slip feet for stability.
Self-lubricating, high-gloss surface helps prevent staining.
Dimensions: 33.97cm x 24.45cm x 2.22cm

(20Z401) Available in all colors



B-Size Safe-Lok® Tray Cover

Fits B-Lok Flat and Divided Trays only.
Side locking mechanism securely locks cover to tray.
Dimensions: 35cm x 25.4cm x 2.775cm

(20Z445)



B-Lok Divided Tray

Holds 12 hand instruments.
Bur holder for 9 FG, RA, or HP burs.
Comes with Insti-Grip.
Non-slip feet for stability.
Self-lubricating, high-gloss surface helps prevent staining.
Holder for Syringe Tip or Saliva Ejector.
Dimensions: 33.97cm x 24.45cm x 2.22cm

(20Z451) Available in all colors



B-Size Safe-Lok® Cover Tall


The Tall size accommodates our Performance Tool Organizer and almost any stainless steel cassettes.
Fits B-Lok Flat and Divided Trays only.
Side locking mechanism securely locks cover to tray.
Dimensions: 35cm x 25.4cm x 3.5cm

(20Z446)



Our Performance Tool Organizer™ fits on Zirc Trays with the Safe-Lok® Tall Cover! See page 13 for E-Z Jett® Cassettes.

FOR B-LOK DIVIDED TRAY



Insti-Grip

To be used with the B-Lok Divided Tray only.
Holds instruments quietly and safely in place.
Dimensions: 16.83cm x .64cm x .95cm

(20Z454I)



Mini Tray

Self-lubricating, high-gloss surface helps prevent staining.
Non-slip feet for stability.
Wide edges for secure grip.
Mini Trays are great for hygiene/exam setup.
Dimensions: 23.81cm x 16.19cm x 2.22cm

(20Z101) Available in all colors



Mini-Lok Tray Cover

Side locking mechanism securely locks cover to tray.
Dimensions: 25.08cm x 16.83cm x 2.54cm

(20Z103)



Mini Tray Cover (Non-Locking)

Dimensions: 23.02cm x 17.15cm x 2.22cm

(20Z103)



B-Size Tray Cover (Non-Locking)

Dimensions: 35.56cm x 25.72cm x .79cm

(20Z441)



Multi-Mod Racks

SAVE ON SPACE & Organize Your Office

Tub, tray, and cassette racks are very effective additions to a dental practice. They save on counter space and several tubs or trays with covers can fit into one rack. These racks have non-slip feet so they are very stable on any counter or cabinet surface. With a chrome-plated finish, these racks can fit any dental practice design.

Multi-Mod Rack 8-Place

Holds 8 Trays or 4 Tubs with covers.
Adjustable to hold Trays, Tubs, and Cassettes.
Sturdy metal construction; chrome-plated.
Non-slip feet.
(21Z410)



Multi-Mod Rack 6-Place


Holds 6 Trays or 3 Tubs with covers.
Adjustable to hold Trays, Tubs, and Cassettes.
Sturdy metal construction; chrome-plated.
Non-slip feet.
(21Z105)



Multi-Mod Rack (8-Place)



Multi-Mod Rack (6-Place)



**Adjustable
to fit your needs!**



*Can adjust to hold
Mini Trays and Cassettes*

WHY USE Zirc's Cassettes?

A low-cost investment with a quick return! Individual instrument handling takes time away from treating patients and increases overall processing time. Zirc's cassettes will prolong your instrument's life and save you time without breaking the bank.

In a 2014 Implementation of Cassettes Case Report, Marie T. Fluent, DDS, writes how the dental office benefited from the transition from no cassettes to using them with every procedure: "Drs. Dale and Thomas were both grateful that they chose to implement ... an Instrument Management System for the organization, processing, and handling of their instruments within their practice. They agreed that cassettes improved all aspects of instrument management... Initially, their opinions were varied; however, at the conclusion of the one-year anniversary, no clinical staff could imagine processing instruments in any other fashion. Cassette usage had become organized, efficient, safe, and routine." Read the article in its entirety at: www.dentaleconomics.com.

- Designed for superior water flow and cleaning - open top, round holes, and radius corners
- Keeps instrument setups together through sterilization and procedure processes
- Protects instruments from extensive damage and prolongs their life
- Non-slip feet for stability and safety
- AUTOCLAVABLE
- Available in 16 colors for color coding by procedure

Sterilization Pouch Size

Cassette	Pouch Size
E-Z Jett® 5-Place	13.34cm x 30.48cm
E-Z Jett® 8-Place	18.1cm x 35.1cm
E-Z Jett® 8-Place Slim	18.1cm x 35.1cm
E-Z Jett® 10-Place	25.4cm x 38.1cm

coming soon!

Pouch-It™

Custom-fit sterilization pouches that will eliminate the guess work of what size pouch to use with each Cassette




Learn more at
www.zirc.com/instrumentcassettes



Keep the
workplace
organized and
confusion to a
minimum!

Provide pleasing presentation chairside.
Minimize stress for patients.


AUTOCLAVABLE Instrument Holders



Syringe Tip Holder

Syringe Tip Holder 5pk

Dimensions: 2.54cm x .95cm x .95cm
(50Z926)



E-Z Jett® Cassette 5-place

Holds 5 instruments.
Great for a basic exam set up!
Syringe Tip Holder included.
Covers stay attached to the cassette and fold back easily.
Fits Statim Sterilizer.
Dimensions: 20.32cm x 9.74cm x 3.02cm
(50Z925) **Available in all colors**



Size comparison with the 8-Place E-Z Jett® Cassette.




E-Z Jett® Cassette 8-place

Holds 8 instruments.
Covers stay attached to the cassette and fold back easily.
Syringe Tip Holder included.
Dimensions: 20.15cm x 11.26cm x 2.85cm
(50Z927) **Available in all colors**




E-Z Jett® Slim 8-place

Holds 8 instruments.
Covers stay attached to the cassette and fold back easily.
Compact size. Fits the Statim.
Stackable in Ultrasonic.
Syringe Tip Holder included.
Dimensions: 18.09cm x 9.84cm x 1.58cm
(50Z911) **Available in all colors**



Instrument Mat

Holds 10 hand instruments.
Reverse side holds 4 surgical instruments.
Non-slip, soft material.
Small size allows easy access to instruments.
Organizes instruments in order of procedure use.
Can be cut to fit Mini-Tray.
Dimensions: 17.5cm x 5.08cm x 0.95cm
(20Z464) **Available in all colors**




Instrument Strap

Mini Mat

Holds 8 instruments.
Organizes instruments at chairside and with sterilization.
Non-slip feet for stability.
Dimensions: 9.84cm x 8.89cm x 1.27cm
(20Z480) **Available in all colors**

Instrument Strap 2pk

Locking strap secures instruments.
Dimensions: 9.20cm x 1.27cm x .79cm
(50Z917) **Gray Only**



E-Z Jett® Cassette 10-place

Holds 10 instruments.
Accessory compartment in the middle for extra storage.
Covers stay attached to the cassette and fold back easily.
Dimensions: 20.16cm x 17.78cm x 2.85cm
(50Z930) **Available in all colors**



Steri-Containers



Closed

Performance Tool Organizer

A deeper, roomier cassette made for large and specialty instruments. Design your own layout with our customizable, removable feet (also sold separately). Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Ideal for:

- Surgical elevators and forceps
- Orthodontic pliers
- Hemostats
- Scissors
- Any instrument that does not fit in your regular cassette



Open

Performance Tool Organizer

Covers stay attached to the cassette and fold back easily.
Dimensions: 20.8cm x 10.8cm x 3.1cm
(50Z940) Available in all colors



Performance Tool Organizer Replacement Feet

Two Strips of 10
(50Z944) Gray only



Holds your large and specialty tools

E-Z Jett® Cassettes fit on Zirc Trays with locking cover and the Performance Tool Organizer® fits our new Safe-Lok® Cover Tall! Check out our Trays & Safe-Lok® Covers on page 10.




Steri-Container

Holds 10-14 instruments.
Keep instrument setups together through sterilization processes.
Dimensions: 20.63cm x 4.76cm x 4.76cm
(50Z900) Available in all colors



Compact Steri-Container

Holds 8-12 instruments.
Easy open latch, just push down.
Statim 2000 holds 7 Compact Steri-Containers.
Dimensions: 18.09cm x 3.81cm x 3.81cm
(50Z905) Available in all colors



Bur Guards



PROTECT YOUR BURS Steri-Bur Guards




Steri-Bur Guards and Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it can also be positioned at three different angles to give you easy access at chairside.

Steri-Bur Guard

- Fits FG, RA, and HP Burs
- Adjustable height for burs of different lengths (12-Hole Tall Bur Guard not adjustable)
- Holds burs in position for maximum cleaning effectiveness
- Drain holes help prevent rusting of burs



12-Hole Bur Guard

Accommodates burs from 12mm-30mm.
Dimensions: 7.3cm x 1cm x 3.5cm
(50Z406) Available in all colors



22-Hole Bur Guard

Accommodates burs from 12mm-30mm.
Dimensions: 12.2cm x 1cm x 3.7cm
(50Z410) Available in all colors



12-Hole Tall Bur Guard

Accommodates burs from 30.9mm-37.1mm.
Dimensions: 7.3cm x .95cm x 4.3cm
Ideal for nose-cone lab burs
(50Z408) Gray only

A Short Bur Adapter is included to stop any shorter burs from falling out.

Universal Bur Adapter 5pk

Can be used with 12- or 22-Hole Bur Guard. Enables long and short burs to be used in the same bur guard without short burs falling out.
(50Z412)






Endo Guards

A unique delivery system to organize, store, and sterilize endodontic files and reamers. The Endo Guard swivels 360° with stops at every 45° for operator's choice. Large slotting provides complete visibility and draining during cleaning.

Steri-Endo Guard

- Holds most lengths for staging, cleaning, and sterilizing files and reamers
- Millimeter gauge measures and positions rubber stopper
- Counter on locking bar tracks the number of sterilizations
- Right- or left-hand use; optional angling
- Available in Hand or Engine styles
- Single-use foam insert cleans/holds files (packaged with 12)



Steri-Endo Guard Engine

Organizes up to 11 engine and 5 hand files/reamers.
Dimensions: 13.65cm x .95cm x 5.56cm
(50Z452) **Available in all colors**

Steri-Endo Guard Hand

Organizes up to 16 hand files/reamers.
Dimensions: 13.65cm x .95cm x 5.56cm
(50Z450) **Available in all colors**



Foam Inserts

Single use.
(50Z451) **48pk**
(50Z451V) **144pk**



Unique organizer rests on your finger.

Endo Assist

Built-in millimeter gauge measures and positions rubber stopper. For right- or left-hand use. Small and large rings are included to accommodate most finger sizes. Flat bottom enables the Endo Assist to stand freely on its own. Single-use foam insert cleans/holds files (packaged with 12).
Dimensions: 3.81cm x 4.45cm x 4.13cm
(50Z460) **Available in all colors**



Assist Stand

Dense foam cleans bio-debris faster and more easily. Holds Endo Instruments in an upright position. Square shape fits great in procedure tub. Thicker foam prevents solution from leaking. Rounded corners for easy cleaning.
Dimensions: 4.76cm x 4.76cm x 3.16cm
(50Z462) **Black Only**



Assist Stand Foam

Single use.
(50Z466) **48pk**
Dimensions: 26.67cm x 19.685cm x 2.22cm
(50Z466V) **144pk**
Dimensions: 26.67cm x 19.685cm x 2.22cm



Bur Blocks



Magnetic Bur Blocks

- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- Non-slip feet for stability
- Drain holes help prevent rusting of burs



8-Hole Bur Block

Covers sold separately.
 Dimensions: 3.81cm x 3.81cm x 1.27cm
 (50Z400) Available in all colors



14-Hole Bur Block

Covers sold separately.
 Dimensions: 7.46cm x 3.81cm x 1.27cm
 (50Z403) Available in all colors



Bur Block Covers

- Covers available to maintain cleanliness of burs during storage
- 1" covers for FG and RA burs; 2" cover (14-hole only) for lab and surgical burs
- Covers are not autoclavable
- Surface disinfect only



8-Hole Bur Block Covers 1" High - Standard

Dimensions: 4.13cm x 4.13cm x 2.54cm
 (50Z401)



14-Hole Bur Block Covers 1" High - Standard

Dimensions: 7.94cm x 4.13cm x 2.54cm
 (50Z404)



14-Hole Bur Block Covers 2" High - HP Lab Burs

Dimensions: 7.94cm x 4.13cm x 5.72cm
 (50Z405)

Rings & Tape


How to effectively identify your instrument setup.

- 1 color for procedure:**
Vibrant Orange is for identifying the procedure.
- 1 color for practitioner:**
Vibrant Blue shows ownership to the practitioner.

Did you know?

The diagonal setup of the orange Rings shows the order of how the instruments should be used. This eliminates confusion between the doctor and assistant.

Place your instrument to determine which ring size you need.



E-Z ID Rings

- Prepackaged in an easy-to-use dispenser
- Available in Small, Large, and XL sizes to fit all instrument handles

25pk

Available in all colors.

- SM (70Z100)
- LG (70Z200)
- XL (70Z000)

Jewel System 200pk

Assorted colors (E, J, T, I, M, and U).

- SM (70Z115)
- LG (70Z215)
- XL (70Z015)

Vibrant System 200pk

Assorted colors (N, O, P, Q, R, and S).

- SM (70Z110)
- LG (70Z210)
- XL (70Z010)

Classic System 200pk

Assorted colors (A, B, D, G, and I).

- SM (70Z105)
- LG (70Z205)
- XL (70Z005)



E-Z ID Tape

- Tape dispenser has built-in tape cutter.
- Dispenser roller prevents tape from twisting.



E-Z ID Tape Roll

Available in all colors.
3m roll.
(70Z300)



E-Z ID Tape System - Jewel

Assorted colors (E, J, T, I, M, and U).
Includes : 0.91m rolls.
(70Z315)

E-Z ID Tape System - Vibrant

Assorted colors (N, O, P, Q, R, and S).
Includes : 0.91m rolls.
(70Z310)

E-Z ID Tape System - Classic

Assorted colors (A, B, D, G, and I).
Includes : 0.91m rolls.
(70Z305)

Label Containers!
Keep containers organized
by labeling them per content.



ISOLATION360™ POWER PACK


Zirc has the most extensive and innovative line of isolation products in the industry. There are practices that may be using one isolation device, not knowing there is something else that could be more beneficial to the procedure. With the Isolation360™ Power Pack, practices are able to explore and utilize different isolation products finding the best fit for their needs. We have an isolation device for every situation!



Isolation360™ Power Pack

Includes:

- 1 Ideal Usage Chart →
- 5 Lingua-Fix®
- 20 SE Cushion
- 10 Pink Petal®
- 9 Airway Armor® (3 Small, 3 Medium, 3 Large)
- 10 Mr. Thirsty® One-Step (5 Pedo/Small & 5 Med/Large)
with one of each HVE adaptor (10mm and 16mm)
- 10 Insti-Dam® (5 Original & 5 Relaxed Fit)
(50Z950X)



See pages 22-27 for individual product descriptions




WHAT PRODUCT SHOULD I USE?

*Product image shows ideal option
○ = Suggested Alternative

	RESTORATIVE				ENDODONTICS		ORAL SURGERY		HYGIENE		ORTHODONTICS							
	Crown & Bridge Inlays/Onlays	Crown & Bridge Inlays/Onlays	Composite / Amalgam	Digital Imaging	Veneer Prep	Veneer Insert	Anterior Tooth	Posterior Tooth	Extractions	Implant	Prophy	Scaling & Root Planning	Sealants	Ultrasonic Scaling	Bonding-Single Bracket	Full Mouth Placement	Removal	Whitening (In-Office)
AIRWAY ARMOR		🦷			🦷	○	○	🦷	🦷					🦷		🦷		
SE CUSHION	○	○								📦						📦		📦
Pink Petal										🦷	🦷	🦷	🦷					
INSTI-DAM			○			📦	📦											
LINGUA-FIX		○ *for cementing		🦷	○							○						
mr.thirsty	🦷	🦷	○		🦷							🦷						

To see the full usage chart, visit...
www.zirc.com/360



A small solution to a BIG problem!

Dental items have been reported as the second most commonly ingested/aspirated foreign objects by adults.

Photo by:
Jesus Escareno





Aspiration & Ingestion Prevention | Patient Protection

There are a multitude of small items used in dentistry: burs, files, implants and veneers to name a few. Combine these items with the moist environment of the oral cavity and you have a realistic probability of something being aspirated or ingested. Airway Armor® was created to eliminate the unnecessary risk for patients and protection for the practice.

Ideal for Crown & Bridge, Endodontics, Extractions, Implants, Inlays/Onlays, Orthodontics, Veneers

Aspiration & Ingestion prevention.
Patient protection & comfort.
Easy placement & removal.
3 convenient sizes - Small, Medium & Large.
Single Use.
Latex Free.



15pk Trail Kit

Includes 5 of each size:
(5) Small
(5) Medium
(5) Large
(50Z960TRI)

Airway Armor® 25pk

Small (50Z960SM)
Medium (50Z960MD)
Large (50Z960LG)

To see an in-use video, visit...
www.zirc.com/airwayarmor

Let your Saliva Ejector Blossom!



Created by hygienist, Karrie Velky, the Pink Petal® is a great isolation product for all dental practices looking to fulfill the need of hands-free isolation. Easily slide the Pink Petal® on any saliva ejector creating a device that will provide continuous suction. This product is great relief for patients who want to hold the suction themselves or choose not to swallow during a procedure.

- Ideal for Scaling & Root Planing, Ultrasonic Scaling, Sealants**
- Holds suction tip in place.
 - Protects & retracts cheek.
 - Patient can close around suction with ease.
 - Easily slides onto saliva ejector.
 - Soft, comfortable material.
 - Single Use.
 - Latex Free.



Pink Petal® 50pk
(50Z975)



To see an in-use video, visit...
www.zirc.com/pinkpetal



The Efficiency of Hands-Free Isolation

Hands-free retraction, isolation and evacuation provided by one simple device! Attach Mr.Thirsty® to your existing HVE, perform the procedure, and dispose when done! Rather than time spent on sterilizing or tracking number of uses, it can be spent on additional procedures, daily tasks or catching a break! According to a time and motion study, Mr.Thirsty® can reduce chair time up to 29%!

- Built-in bite block for patient comfort
- Connects directly to HVE Valve
- Keeps work environment dry
- Fully Assembled
- Latex-Free



WITH ONE STEP, ACHIEVE
SIMPLICITY

AFFORDABLE

no upfront investment

DISPOSABLE

no time wasted in maintenance or sterilization

ADAPTABLE

connects to existing HVE, allowing it to be used in any operator

COMFORTABLE

provides tongue and cheek retraction and includes a built-in bite block, creating comfort for you and your patient

TRIMMABLE

available in two sizes that can be trimmed to fit, no need to overstock with a variety of sizes

mr. thirsty® one-step



REDUCE CHAIR TIME
BY AS MUCH AS

29%



Take comfort to
the next level with
lightweight tubing!



10pk Trial Kit

10mm Adaptor (50Z98710)
16mm Adaptor (50Z98716)
Includes 5 of each size:

- (5) Pedo/Small Adult
- (5) Medium/Large



25pk

- Pedo/Small Adult
10mm Adaptor (50Z986SM10)
16mm Adaptor (50Z986SM16)
- Medium/Large
10mm Adaptor (50Z986LG10)
16mm Adaptor (50Z986LG16)

ALSO
AVAILABLE IN
100 PACK

100pk

- Pedo/Small Adult
10mm Adaptor (50Z986SM-B10)
16mm Adaptor (50Z986SM-B16)
- Medium/Large
10mm Adaptor (50Z986LG-B10)
16mm Adaptor (50Z986LG-B16)



Comfort Kit

The Comfort Kit's lightweight tubing allows your heavy HVE Valve/Tubing to stay mounted to the equipment.

Includes:

- Rotary plastic HVE Valve
- 1.83m lightweight tubing
- mr. thirsty® one-step NOT included (50Z992)



10mm/16mm Adaptors

10mm Adaptor 1pk (50Z993)
16mm Adaptor 1pk (50Z994)

136
SSS



Dental PRODUCT SHOPPER
RECOMMENDED
2014


To read a full product evaluation visit...
<http://www.dentalproductshopper.com/relaxed-fit-insti-dam>

Easy to use, non-threatening,
and comfortable for patient.

INSTI-DAM®

Features for BOTH Relaxed Fit and Original Insti-Dam®

- Built-in flexible medium gauge material that comes pre-assembled with frame and pre-punched hole
- Pre-punched hole makes positioning easy and allows eyes and nose to be open and comfortable
- Helps eliminate tearing (additional holes can be punched)
- Compact design fits outside patients' lips, easy-to-use, non-threatening and comfortable for patients
- Radiographs may be taken by folding it to the side
- Single use



Latex-Free Insti-Dam® Relaxed Fit 20pk

Ideal for posterior isolation.

The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing, or pulling.

Dimensions: 10.79cm diameter

Patent #7044738

(50Z457)

Blue Latex-Free Insti-Dam® 20pk

Ideal for anterior isolation.

Made with polyisoprene.

Dimensions: 10.79cm diameter

(50Z459)

Natural Latex Insti-Dam® 20pk

Ideal for anterior isolation.

Dimensions: 10.79cm diameter

(50Z455)

Lingua-Fix® & SE Cushion



Lingua-Fix®

Lingua-Fix® provides a dry work area.
Holds tongue in a steady position to protect from injuries.
Simple operation...connects to low volume suction.
Comfortable, no sharp edges or corners.
Large suction surface avoids blockage and is more comfortable.
Reduces sedimentation of drilling by-products in pipes and tubing.
Single use .
Comes pre-assembled.
(50Z953S) **10pk**
(50Z953) **50pk**



SE Cushion

Used in conjunction with the Saliva Ejector.
Prevents tissue grabbing and tissue irritation.
Increases patient comfort.
Improves range of suction.
Quiets the noise of the suction.
Single use.
(50Z970) **100pk**



HELPFUL TIP:
Ideal for use
with rubber
dams

**For comfortable and efficient
suction or evacuation.**

The Relax Factor:

The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing or pulling.



Insti-Dam® Dispenser

Made for Insti-Dam® Original ONLY
Can be mounted to wall or cabinet or can sit on a counter.
Available in White (A) for latex or Neon Blue (N) for latex-free.
Non-slip rubber bottom for stability.
Holds 35 Insti-Dams®.
Dimensions: 12.38cm x 12.38cm x 15.56cm
(50Z471) **Color Code: A & N**

AM S



High Definition Imagery




CRYSTALHD®

Experience high-definition quality. The brighter the workspace, the better the illumination and visibility. With ZIRC's Crystal HD® lens, you are visually working in a 40% brighter environment compared to the Rhodium lens that is currently the industry standard. If you haven't made the change to HD - what are you waiting for?

**40%
BRIGHTER**
Than Rhodium
Mirrors



CUTTING EDGE IMAGERY



Rhodium = Reflects color 75% accurately
Crystal HD® = Reflects color 99% accurately

Provides exceptionally bright,
shadow-free images!



Crystal HD®



Available in multiple styles!

Designed to reduce eye strain and fatigue.

Autoclavable



As a microscope user, I can understand the importance of clear vision. ZIRC mirrors significantly improved the brightness of my image, giving my eyes highest quality vision. It's true - brightness makes all the difference.

- Dr. Dimitris Tsanaktsidis, Microscopic & Aesthetic Dentistry

CRYSTAL HD® MIRROR + MIRROR MAGIC®
The Perfect Pair

Pair our Crystal HD® lens with our mirror magic® System for the highest quality vision and best procedural outcome. See page 32.




Crystal HD® Mouth Mirrors & Handles

Flush surface, no debris build-up can occur



Resin


Gap where debris build-up occurs



Stainless Steel

	RESIN	STAINLESS
Autoclavable	X	X
Antimicrobial protection	X	
Lens flush with mirror frame	X	
No galvanic shock on patient	X	
Rust-free	X	

Why use Zirc's Resin Mirrors?



Assorted Soft Grip 12pk

Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

#4 (50Z363-NEON)

#5 (50Z364-NEON)

Jewel Assorted

2 of each color: E, J, T, I, M, and U.

#4 (50Z363-JWL)

#5 (50Z364-JWL)



Assorted Thin Grip 12pk

Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

#4 (50Z353-NEON)

#5 (50Z354-NEON)

Jewel Assorted

2 of each color: E, J, T, I, M, and U.

#4 (50Z353-JWL)

#5 (50Z354-JWL)



Black/Gray Soft Grip 12pk

#4 (50Z363)

#5 (50Z364)

Black/Gray Soft Grip 4pk

#4 (50Z361)

#5 (50Z362)



Handles




Ergo Grip Handle (Simple Stem)

(50Z347) Color Code: B-S, and K (black)



What's your mirror size?



#0 - 14mm

#3 - 20mm

#4 - 22mm

#5 - 24mm

diagram is to scale



Resin Screw-on Heads 12pk

Flush surface, no debris build-up can occur.
Simple stem.

1-Color 12pk

Available in all colors and K (black).

#4 (50Z392)

#5 (50Z393)

Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

#4 (50Z392-NEON)

#5 (50Z393-NEON)

Jewel Assorted

2 of each color: E, J, T, I, M, and U.

#4 (50Z392-JWL)

#5 (50Z393-JWL)

Classic Assorted

A, B, D, G, I, and K (black).

#4 (50Z392-PSTL)

#5 (50Z393-PSTL)



Vibrant Assorted Double-Sided Mirror Head 4pk

1 of each color: N, O, R, and S.
Simple stem.

#4 (50Z394-NEON)



Gray Double-Sided Mirror Head 4pk

Simple stem.

#4 (50Z394)



Stainless Steel Head 12pk

Simple stem.

#3 (50R389)

#4 (50R390)

#5 (50R391)



Stainless Steel Head 1pk

Simple stem.

#0 (50R388)



Zirc's Mirrors & Handles are autoclavable!

mirror magic® Anti-Fog System



without mirror magic®



with mirror magic®



mirror magic®
anti-fog system



Increase Your Chairside Efficiency

The mirror magic® Anti-Fog System prevents fogging on mouth mirrors and offers the practitioner a quick method for removing debris during procedures. The disposable Swipe Pads stick to the back of a gloved hand or a patient bib for easy access. The Anti-Fog System requires only a small amount of solution and is economic in size, making it easy to store.

mirror magic® Trial Kit

2-6ml bottles of Anti-Fog Solution
24 Swipe Pads
(50R200X)



Watch the magic in action
at www.zirc.com/mirrormagic

Combine mirror magic® with Crystal HD®
Mouth Mirror for the best view in dentistry.





"The mirror magic® Anti-Fog System decreases fogging while also decreasing the build-up of debris on my mirror. This affords me outstanding visualization, which is crucial to performing quality dentistry. mirror magic® has truly made a difference in my practice and is one of those features I didn't realize I needed; and now I can't imagine practicing without it."

— Andrew M. Goldsmith, DDS, DICOI, FIALD



10pk Anti-Fog Solution (6ml ea.)

1 bottle lasts approx. 50 patients.
(50R205)



Swipe Pads 96pk

Single use.
(50R210)



mirror magic® Organizer

Holds: 48 Swipe Pads and 1 - 6ml bottle of Anti-Fog Solution.
Store on operatory counter or in drawer.
Non-slip bottom for secure placement.
Note: box only, does not include mirror magic.
Dimensions: 3-3/4" x 3-5/8" x 1-7/8"
(50R201)

S



\$122 Value!



Clarity Collection

Includes:
Anti-Fog Solution 10pk
Swipe Pads 96pk
1 Crystal HD® Mouth Mirror
1 mirror magic® Organizer
(50R220)

The Perfect Pair

The superior visual acuity and 40% brighter Crystal HD® lens combined with the quick debris removal and Anti-Fog Solution of mirror magic® have proven to be the perfect pair. Increase chairside efficiency and speed up your process with this dynamic duo.





CAD CAM BLOCK LOCKER™

Features

- Stores and manages CAD/CAM milling blocks and accessories
- Drawer accommodates accessories such as burs & lubricants
- Accommodates all blocks
- Wall mounted unit features a dark gray powder-coated finish for a modern look
- Provides a professional platform to present and discuss same-day restorations
- Online customizable tags to create organization by color and shade (<https://www.zirc.com/products/cad-cam-block-locker>)



CAD/CAM Block Locker™ Large

Keep blocks in original packaging.
 Holds up to 440 blocks in a combination of packages and individual blocks.
 Dimensions: 63.5cm x 67.63cm x 23.49cm
 Weight: 22.41kg
 (50Z416) Color Code: E, J, N, Q, S, T, U

CONTROL BLOCK CHAOS
 SHOWCASE YOUR MILLING BLOCKS IN STYLE



CAD/CAM Block Locker™ Small

Designed to store individual/loose blocks only.
 Holds up to 288 individual blocks.
 Dimensions: 49.53cm x 40.64cm x 23.49cm
 Weight: 13.03kg
 (50Z415) Color Code: E, J, N, Q, S, T, U

Pick Your Color!



Denture/Retainer Boxes



Retainer Box

One piece construction.
Solid locking mechanism.

Retainer Box 2.54cm Deep 12 pk

Dimensions: 7.62cm long x 2.54cm deep

(25X500) 12pk

Color Code: J, S, R, N, O, P, Assorted

Retainer Box 3.81cm Deep 12 pk

Dimensions: 7.62cm long x 3.81cm deep

(25X550) 12pk

Color Code: J, S, R, N, O, P, Assorted

Solid locking mechanism so your retainer stays put!

S

Great colors to choose from!

Denture Box

One piece construction.
Solid locking mechanism.

Denture Box 12 pk

Dimensions: 9.53cm long x 4.45cm deep

(30X800) 12pk

Color Code: T, U, E, I, J, M, Assorted

S



Saliva Ejector Accessories & MEGA Products



Saliva Ejector Valve (with Shut-Off)



Quick disconnect with swivel.
Rubber tip with protective screen.
Provides variable control from full on to full off.

Saliva Ejector Valve with Lever On/Off Control

Dimensions: 7.62cm x 1.59cm x 1.91cm

(23E365) Gray Only

Saliva Ejector Valve Lever Only

Dimensions: 1.59cm x .48cm x 1.59cm

(23E364) Gray Only

Saliva Ejector Replacement Tip

Dimensions: 1.59cm x 1.59cm x 2.54cm

(23E363) Gray Only

Hands free isolation!
See page 20-27.



Saliva Ejector Screens

Single use.

Saliva Ejector Screens

Dimensions: .95cm Diam. x .95cm

(50Z952S) 25pk

(50Z952P) 250pk



MEGA Products

MEGA Gracey

39.37cm in length.
Great teaching tool.
A gift for any dental personnel.
(50Z312) Mustard Only



MEGA Mirror

(Patient hand mirror).
6.35cm lens and 39.37cm in length.
Keeps children entertained.
Great gift for any dental personnel.
(50Z309) Color Code: N, P, Q, R, S



Fun gift ideas for any
dental personnel.
Keeps children
entertained.



Water Wise & Steri-Soaker



Water Wise Starter Kit

Holds 1 liter.
Includes bottle, adaptor, lid, neon yellow and neon blue bubbles (water level indicator), and tube.
(23Z100)



Water Wise Replacement Bottle

Holds 1 liter.
Includes bottle, blue lid, and neon yellow and neon blue bubbles (water level indicator).
(23Z105)



Water Wise vs. Traditional Water Bottle

	Water Wise	Water Bottle
Easy to clean	X	
Antimicrobial	X	
Easy-to-fill large mouth	X	
Quick release	X	
Level indicator color coded for: Sterile and distilled water, Medicament	X	
Strong	X	X
Easy to grip	X	
Sealable lid	X	
No threads to strip	X	
Pressure relief valve	X	

We strive to make things easier around the dental office.

The Steri-Soaker

An efficient, high-volume germicide tray designed for maximum use of chemicals to clean/disinfect instruments. It is designed for easy cleaning and is large enough to accommodate ZIRC Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.

Steri-Soaker

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters.

Hinged fitted cover.

Safe-to-use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines.

Dimensions: 35.24cm x 24.45cm x 13.02cm

(50Z805) Color Code: B, I



Designed for easy cleaning!



E-Z Access Shelf/Accessories

E-Z Access Shelf

E-Z mounting to standard equipment post (4.45cm-6.98cm); No tools needed.
Will not damage post's finish.

E-Z adjustment to multiple positions.

Contoured edge to contain spills and materials.

Stain resistant finish; surface disinfect only.

Weight capacity of 11.33kg

Workable tray area: 34.93cm L x 26.04cm W.

Standard E-Z Access Shelf as is accommodates 4.45cm-5.72cm post.

Rotates 360°.

(20Z420) Color Code: A, G, I



Provides convenient work space for auxiliary items.



E-Z Access Shelf Adapters

E-Z Access Shelf Adapter for Smaller Posts

Adapts to posts 3.49cm-6.99cm

Dimensions: 5.13cm x 4.76cm x 1.91cm

(20Z422) Color Code: A, G, I

3-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 5.72cm-6.03cm

Dimensions: 10.16cm x 5.08cm x 5.08cm

(20Z437) Color Code: A, G, I

4-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 6.03cm-6.99cm

Dimensions: 11.43cm x 5.08cm x 5.08cm

(20Z438) Color Code: A, G, I



Cup Holder-Snap On

Fits all standard cups.

No tools needed.

Dimensions: 7.94cm x 6.35cm x 2.54cm

(20Z436I)



E-Z Access Shelf Surface Mount

Mount on horizontal flat surface only.

Dimensions: 9.53cm x 9.53cm x 7.30cm

(20Z424) Color Code: A, G, I



Our products are able to fit every need!



8" Extension Arm

Enables shelf to extend over patient.

Shelf has strength to hold 25 lbs. when connected to arm.

May be used with Wall Mount, Surface Mount, or directly to post.

Dimensions: 25.4cm x 10.16cm x 11.43cm

(20Z428) Color Code: A, G, I



Wall or Cabinet Mount

Allows you to mount the E-Z Access Shelf onto a wall or the side of a cabinet. 90° angle of mount permits 180° movement of shelf on cabinet or next to wall.

Dimensions: 16.51cm x 14.61cm x 12.38cm

(20Z426) Color Code: A, G, I



Cabinet Trays & Lab Pans

AM S 

Cabinet Trays



No 16A – Hand Instruments (Regular Well)

Dimensions: 20cm x 9.53cm x 2.38cm
(20Z201) Color Code: A, B




No 16 – Hand Instruments (Deep Well)

Dimensions: 19.69cm x 9.53cm x 2.54cm
(20Z202) Color Code: A, B



No 17 – Clamps, Retainers, etc.

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z203) Color Code: A, B



No 18 – Finishing Strips

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z204) Color Code: A, B



No 19 – Pliers, Scissors, etc.

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z205) Color Code: A, B



No 20 – Discs, Wheels, Stones

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z206) Color Code: A, B



Wall-Hanging



AM S 

Lab Pans

Organizes and transports impressions, models, dentures, and crown and bridge appliances.
Stainless-steel spring clip holder.
Color-coded.

Lab Pan

Dimensions: 19.05cm x 14.61cm x 17.78cm
(17Z101) Color Code: A, B, D, G, M

Wall-Hanging Lab Pan

Same as standard version but with strip on top back for hanging on any vertical surface.

Saves on installing expensive shelving.
Dimensions: 19.69cm x 14.61cm x 17.78cm

(17Z102) Color Code: A, B, D, G, M



Wall Strip Lab Pan

Drilled and countersunk with screws and wall anchors.
45.72cm Wall Strip.

Holds 3 Lab Pans.

Dimensions: 45.72cm

(17Z104)

S

Save on installing expensive shelving.



Deluxe Cotton Roll Dispenser

Heavy-duty no-skid base.
Dispenses up to three rolls at a time.
Extra sturdy construction.
Dimensions: 10.64cm x 5.24cm x 10.16cm
(28R800) **Color Code: A, B**

S



Cotton Roll Holder

Stainless pins embedded in weighted bottom.
Neatly holds cotton rolls in an upright position.
Hinged lid helps keep cotton rolls clean, prevents spilling.
Dimensions: 9.21cm x 8.26cm x 6.03cm
(28R900) **Color Code: A, B**

S



Pellet Dispenser

Spring loaded – pushes pellets to the surface.

Dispenser for Small Pellets (sizes 3, 4, 5)

Dimensions: 4.76cm x 4.76cm x 5.24cm
(28R825) **White Only**

Dispenser for Medium Pellets (sizes 0, 1, 2)

Dimensions: 4.76cm x 4.76cm x 5.24cm
(28R850) **White Only**

S

Sponge Dispensers

Spring action releases pads one at a time.
Bottoms are weighted to prevent tipping.

5.08cm x 5.08cm Sponge Dispenser

Dimensions: 6.35cm x 6.35cm x 7.3cm
(28R950) **Color Code: A, B**

10.16cm x 10.16cm Sponge Dispenser

Dimensions: 11.43cm x 11.43cm x 7.3cm
(28R975) **Color Code: A, B**

S



Matrix Dispenser

One roll of 1cm x 11m Matrix clear tape 0.002 gauge.
Allows Matrix to be contoured (u-shaped) as it dispenses.
Stainless steel cutter allows for desired length.
Disposable.

Dimensions: 5.08cm x 5.39cm
(55Z200) **Color Code: A, B**

S

E-Z Storage/Handi-Hopper

E-Z Storage Tub System

Hinges allow cover to stay open on its own!
Each Tub comes with one vertical divider and one half-depth adjustment divider.
Hinged cover comes in white or clear.

E-Z Storage Organizer (with white cover)

Dimensions: 14.92cm x 12.38cm x 14.61cm
Rack sold separately.
(29R605)

E-Z Storage Organizer (with clear cover)

Dimensions: 14.92cm x 12.38cm x 14.61cm
Rack sold separately.
(29R610)

Replacement Hinges (pair)

Dimensions: 2.22cm x 1.27cm x .32cm
(29R627)

S



E-Z Storage Single Rack

Holds storage organizer Tub.
Tub can be placed flat or at an angle for easy access of supplies.
Wall or cabinet mounted with optional mounting bracket.
Sturdy metal construction; plastic-coated.
White, non-slip feet.
Dimensions: 17.78cm x 12.7cm x 12.06cm
(29R650)

S

E-Z Storage Mounting Bracket

White metal bracket for the E-Z Storage Tub single rack.
Pre-drilled holes for cabinet or wall mount.
Dimensions: 17.78cm x 12.7cm x 12.07cm
(29R651)

S



AM S

Mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste.

Handi-Hopper (with Mounting Bracket)

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean.
Slotted mounting bracket (included) for easy removal and replacement.
To use for infectious waste, order red (Biohazard) liners.
To use for other waste, order clear liners.

Handi-Hopper

Dimensions: 17.78cm x 10.16cm x 16.19cm
(20Z404) Color Code: G, I

Handi-Hopper Bracket

Dimensions: 5.72cm x .64cm x 2.54cm
(20Z408) Color Code: G, I



Handi-Hopper Liners (Single-Use)

Handi-Hopper Liners 100/Box (Red Biohazard)
(20Z405) 100pk

Handi-Hopper Liners 100/Box (Clear)
(20Z410) 100pk

AM S



Waste-Drop Ring

Waste-Drop Ring for cabinet top.
Fits 12.7cm hole.
Made of self-lubricating, easy-to-clean resin.
Dimensions: 14.29cm x 14.29cm x 7.62cm
(22Z100) Color Code: A, G, I

Cuspidor/Dispos-A-Screen Traps & Needle Capper



Reduce possibilities of infection and contamination.

Cuspidor Screens AM

Single use.

Cuspidor Screens Large Universal – 6.35cm (gray)

Adec, Ampco, Belmont, Dab Atlante, Dansereau, DCI Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter, Schein.

(51Z205S) **Gray 25pk**

(51Z205V) **Gray 144pk**

Cuspidor Screens Small – 4.13cm (white)

Adec, Cascade and Performer cuspidor bowls.

(51Z200S) **White 25pk**

(51Z200V) **White 144pk**

Dispos-A-Screen Traps

Avoid contact with potentially hazardous material. Improves performance of suction system. Color coded by size. Single use.

(25pk)

(144pk)

Blue – 4.5cm

Adec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus MDT, Schein.

25pk (51Z045S), **144pk** (51Z045V)

Green – 4.76cm

Ampco, Belmont, Biotec, Chayes, Dabi-Alante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma and others.


25pk (51Z123S), **144pk** (51Z123V)

Gray – 7cm

Midwest 210 Models.

25pk (51Z070S), **144pk** (51Z070V)

Left to right: Blue, Green, Gray, Yellow, Beige, Mauve and Coral



Yellow – 5.72cm

Adec, model with metal canisters.

25pk (51Z040S), **144pk** (51Z040V)

Beige – 7.3cm

Den-Tal-EZ models with metal canister.

25pk (51Z060S), **144pk** (51Z060V)

Mauve – 3.18cm

Adec, Cascade, and Performer.

25pk (51Z080S), **144pk** (51Z080V)

Coral – 6.67cm

Pelton and Crane Spirit Unit with Plastic Canister.

25pk (51Z050S), **144pk** (51Z050V)



Needle Capper

One-handed capping and recapping.

Syringe is for use; at convenient position in Needle Capper.

Base of device is secured to tabletop or tray with adhesive strip (included).

Needle is uncovered at same time syringe is withdrawn from Capper.

After use, syringe is replaced for recapping.

Needle Capper available in Standard and Long sizes.

Either size tube will fit standard base.

Standard Needle Capper

(base, one tube, adhesive tape)

Dimensions: 5.56cm x 3.18cm x 3.18cm

(27R700) **White Only**

Long Needle Capper

(base, one tube, adhesive tape)

Dimensions: 6.35cm x 3.18cm x 4.45cm

(27R704) **White Only**



Mighty Mixer Bowls, Spatulas & Organizers

Single-use
disposable mixing
bowls.



Mighty Mixer Bowls 36pk (White Only)

Disposable, eliminating cleanup time.
Fits securely into mixing bowl without slippage.
Saves wear on mixing bowl.
Free-standing for use with or without mixing bowl.
Capable for use with electric Alginate mixers.
Dimensions: 12.38cm Diam. x 7.93cm
Single - Use.



Mighty Mixer Bowls 36pk (Assorted)

Same features as the white Mighty Mixer Bowls.
Dimensions: 12.38cm Diam. x 7.93cm
9 of each color: J, N, P, and R
Single - Use.
(50Z505-ASTD) 36pk



Mighty Mixer Spatula 4pk (Assorted)

Same features as the white Mighty Mixer Spatula.
Dimensions: 18.57cm .48cm x 3.18cm
1 of each color: J, N, P, and R
(50Z501-ASTD) 4pk



Mixing Sticks 100pk

Single use.
Dimensions: 11.43cm long
(50Z521) 100pk



Cement Spatula

White plastic, double-end spatula.
For fillings, cement, and pastes.
Self-lubricating plastic for easy cleaning.
Dimensions: 16.19cm x .64cm x .95cm
(50Z502)



Delrin Spatula 10pk

Flexible flat blade.
Non-stick.
Resistant to chemicals.
Great for mixing cements, acrylics, etc.
Dimensions: 11.74cm x .47cm Diam.
(50Z523) 10pk



Mighty Mixer Spatula

Self-lubricating for easy cleaning.
White plastic spatula.
Broad, medium flex blade.
Ideal for alginates and plaster.
Dimensions: 18.57cm .48cm x 3.18cm
(50Z501) White Only



Implant Organizer

Quick and easy recognition. Individual compartments for implants.
Shape of the organizer matches the arch of the mouth for implant placement.
Separate compartments prevent implant appliances from mixing together.
Lubricated plastic enables tooth number to be written in the designated area.
Cover helps keep materials secure.
Dimensions: 17.46cm x 11.75cm x 1.59cm
(20Z490B) Blue Only



No more
confusion!



Restorative Organizer

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays.
Disposable Mixing Wells included for bonding agent and cement.
Accessory area for brushes and materials.
Lubricated plastic enables tooth number to be written in the designated area.
Light protective cover can be used as a light shield and to prevent materials from setting early.
Includes 25pk of Mixing Wells.
Dimensions: 17.46cm x 11.75cm x 1.59cm
(20Z495R) Neon Purple Only



Mixing Well (Single-Use) 144pk (20Z496V) 144pk



MAKING DENTAL TEAMS
HAPPIER & MORE *efficient!*



Because Time Is Everything

www.zirc.com | +1 763.628.6636 | sales@zirc.com

International Product Guide