



CELEBRATING
50
YEARS

Zirc

Because Time Is Everything

PRODUCT GUIDE 2017

www.zirc.com | +1 763.682.6636



WHY ZIRC

2

COLOR CODING

4 - 21

Why Color Coding?	4
Color Chart	5
Tub Kits	7
Tub Accessories	8
Trays & Covers	10
Multi-Rod Racks	12
Cassettes	13
Steri-Containers/Mats	15
Bur Guards/Blocks	16
E-Z ID Rings & Tape	19
Zobe	20

ISOLATION

22 - 25

mr.thirsty® one-step	22
Insti-Dam® SE Cushion Lingua-Fix®	24

MOUTH MIRRORS

26 - 31

Crystal HD® Mirrors	26
mirror magic® Anti-Fog Solution	30

OPERATORY/MISC

32 - 40



Celebrating 50 Years!

This year, we're celebrating 50 years of offering our customers the highest quality products at the best value. James Campion, our CEO, founded Zirc in 1967 and his daughter Nicolle, who was recently promoted to President, joined the team in 2000. Dental technology and processes have seen many changes over the last 50 years, but at Zirc, our commitment to providing time-saving, innovative dental products continues to be our top priority.

Zirc Perks

Members receive FREE samples, money-saving promotions, tips to get organized and stay organized, product evaluation invites, and more! Sign up to become a Zirc Perks Member via www.zirc.com or call us at +1 763-682-6636 (available in U.S.A. and Canada only).



Made in the USA

Since 1967, Zirc has provided the dental industry with new and innovative products used for sterilization and organization in the dental office. Our headquarters, located in Buffalo, Minnesota, houses our engineering, tooling, manufacturing, distributing, and customer service operations. We're proud to say that our state-of-the-art, high-quality products are made in the U.S.A.



Antimicrobial

Antimicrobial compounds are intended to protect parts against staining, odors, and degradation due to fungi, mold, yeasts, mildew, algae, and other microbes. Zirc incorporates antimicrobial compounds during the manufacturing process. Unlike a coating, this prevents fading or wearing away, giving you long-lasting protection. The incorporation of antimicrobial compounds in our plastic products is not intended to protect users or others against food-borne or disease-causing bacteria. For more information, please reference EPA Registration Notice 2000-1.



Experience Our Products* for FREE

The ZOBE kit contains over 25 of our top-selling products for a stress-free, fun, organized practice. ZOBE ships at no cost, allowing you to explore the products before you buy. No charge shipping and return. Only available in USA, contact Zirc for more information. See page 20 for more details.



*for demo/viewing only




WANT TO MAXIMIZE SPACE?

E-Z Jett® Compact 8-place Slim

Our autoclavable E-Z Jett® Cassettes extend the life of your instruments by protecting and keeping them together through your procedures and sterilization process. They are designed for superior water flow during cleaning and have non-slip feet for stability and safety.

- Holds 8 instruments
- Covers attached to the cassette open and fold back easily
- Compact size. Fits the Statim
- Maximizes the space of the Ultrasonic and Autoclave
- Stackable in Ultrasonic. Cassettes lightly attach to one another
- Antimicrobial incorporated during manufacturing process
- Syringe Tip Holder included
- Autoclavable



Size comparison with the 8-Place E-Z Jett® Cassette.

**Maximize space!
Stackable in Ultrasonic
and Autoclave!**

AVAILABLE SPRING 2017

Performance Tool Organizer

A deeper, roomier cassette made for specialty instruments that will not fit in standard cassettes. Design your own layout with our customizable, removable feet. The feet can be positioned to support and protect your valuable instruments. Use the Performance Tool Kit in conjunction with any of our E-Z Jett® Cassettes, which hold your hand instruments. Designed for superior water flow for superb cleaning. See page 15 for full details!

**Holds your
large and specialty
tools**



**Design your own layout
with our customizable,
removable feet**



ORGANIZE WITH EASE Why Color Coding?

Increase Profitability

Color coding decreases time for setup and clean up, allowing you to see more patients. It also reduces the amount of inventory needed, limiting expired materials.

Other Benefits


- Simplified training of new staff
- Decreases amount of inventory – fill one tub vs. stocking multiple treatment rooms
- Easy and efficient re-stocking of dental materials
- Reduces the amount of expired materials
- Reduces the risk of an exposure incident
- Eliminates the need to stop a procedure so the assistant can find a missing instrument or material
- Contains dirty instruments in covered, puncture-proof, leak-proof containers satisfying the CDC and OSHA standards
- Small items stay organized and at your fingertips
- Set-ups appear organized, sanitary, and visually pleasing to patients

Want help on how to color code in your dental office?

For more information, visit us online to download our Color Coding brochure.




Color coding has been proven to save time for each procedure and will allow you to see more patients in a day.



ZIRC Color Chart

Classic




Jewel



NEW COLORS

Vibrant





Syringe Composite Tub



Capsule Composite Tub



Endo Tub



Crown Prep Tub



Hygiene Tub



Implant Tub

SMALL CHANGES Big Results

Cluttered drawers = stress and chaos. Stay efficient with color-coded organization.

Start by removing everything from your operatory drawers. Place in a large box and set aside. Restock your operatory drawers with staple supplies such as: bibs, saliva ejectors, cotton products, etc. (basically, anything you use in EVERY procedure).

Here's where the Tubs come into play. You need a minimum of one Tub for each procedure. Each Tub should be a different color so that you can quickly identify procedures by color. Stock each Tub with all materials you will need for that procedure.

Store ALL Tubs in a central area (such as sterilization). Take the Tub into the treatment room when a procedure is performed. Pull out what you need at chairside. Disinfect and re-store.

**Reduce your stress and have fun
with color organization.**





MAKE IT SIMPLE WITH AN All-In-One Kit

Our kits are designed to store and organize materials and keep your inventory under control. In addition, they are easily portable for transportation between procedures and your centralized location for sterilization, disinfecting, and restocking.


Complete Tub

Includes: Procedure Tub, Tub Liner, Tub Dividers, Cup/Cover Single, Cup/Cover Double, Divided Slide Tray, Safe-Lok Cover for Tub®.

Dimensions: 32cm x 30cm x 10.6cm

(20Z455) Available in all colors

AM S



Use: Any Procedure

Composite Kits

Store and organize composite materials with our convenient Composite Kits. When stored in a Tub, contents become portable and can be transported to a central sterilization area for disinfecting and restocking. Organizational tools help you keep better inventory control with Clip-On Holders that allow you to separate composite for use during procedures. Included Long Tub Cup and Slide Tray provide additional organizational space and Shade Label Sheet provides self-stick labels preprinted with composite shade.

Syringe Composite Kit

Syringe holder allows you to store up to 40 syringe tubes. Kit includes: Syringe 20-Unit Stand, 1 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Syringe Clip-On Holder, Shade Label Sheet 2pk.

(20Z489)

AM S



Capsule Composite Kit

Long Tub Cups provide a total of 20 compartments for capsule storage. Includes: Tub Liner, 4 Long Tub Cup/Cover, 1 Divided Slide Tray, 1 Capsule Clip-On Holder, Shade Label Sheet 2pk.

(20Z479)

AM S



Build Your Own Tub



Procedure Tub

Allows you to “build” your own setup.
 Dimensions: 31.12cm x 27.62cm x 6.98cm
 (20Z463) Available in all colors



Safe-Lok® Tub Cover

Locking cover for securing materials.
 Keep materials protected during storage and transportation.
 Dimensions: 32.39cm x 30cm x 3.5cm
 (20Z459)



Tub Liner

Foam liner that sits in bottom of Tub.
 Tub not included. Replace as needed.
 Dimensions: 24.45cm x 22.54cm x .32cm
 (20Z462)



Tub Box

Designed to hold/organize “hard-to-find” supplies.
 Non-slip bottom for secure placement.
 Dimensions: 21.27cm x 5.87cm x 2.86cm
 (20Z474) Gray Only



Tub Insert

10 compartments for easy organization.
 Sits on rim of Tub, allowing for storage underneath.
 Dimensions: 27.3cm x 24.6cm x 2.54cm
 (20Z465) Gray Only



Tub Dividers

Creates up to 24 compartments.
 May be cut or trimmed to create larger compartments.
 Dimensions: 32.39cm x 30cm x 3.5cm
 (20Z461) Gray Only



Use color-coded tape on Tub Cups and Tub Box to identify procedural contents. (Learn more on page 19)



Tub Cups with Covers

Fill Tub Cups with brushes, clamps, cotton pellets, prophylaxis angle tips, prophylaxis paste, pumice, syringe tips, implant tools, etc. Can be used in procedure tubs, drawers, on trays, etc. for small, easy-to-lose items. Covers are transparent.



Single Tub Cup with Cover

Dimensions: 5.27cm x 4.13cm x 2.86cm
 (20Z471) Gray Only

Double Tub Cup with Cover

Dimensions: 8.6cm x 5.72cm x 2.86cm
 (20Z472) Gray Only

Long Tub Cup with Cover

Dimensions: 21.27cm x 5.72cm x 2.86cm
 (20Z473) Gray Only



Divided Slide Tray

Dimensions: 24.77cm x 13.65cm x 2.7cm
(20Z466) Clear

Flat Slide Tray

Dimensions: 24.77cm x 13.65cm x 2.7cm
(20Z467) Clear



Shade Label Sheets 2pk


Preprinted labels for marking Tub materials by procedure. Shade Labels can be used with the Syringe Composite and Capsule Composite Tub Organizers.
(20Z463-LBL)

Slide Trays

Sits on rim of Tub and allows for storage underneath.



Check out our Implant and Restorative Organizers on page 40.



Syringe 20-Unit Stand

Stands up for easy access to composite.
Lays flat to fit inside Tub or drawer.
Holds up to 20 syringe tubes.
Has space for 2 Syringe 10-Unit Add-Ons (Stand can hold up to 40 syringes).
Dimensions: 22.22cm x 12.06cm x 3.96cm
(20Z487)



Syringe 10-Unit Add-on

An Add-on for the Syringe 20-Unit Stand.
Gives additional space for syringe storage.
Holds 10 syringe tubes.
Dimensions: 12.74cm x 5.87cm x 1.9cm
(20Z488)



The Syringe 10-Unit Add-on may be added for storage of up to 40 syringe tubes.



Syringe Clip-On Holder

Allows quick access to composite during procedure.
Conveniently clips on any side of Tub, B-Lok Tray or Capsule Composite.
Holds 15 syringe tubes and 3 applicator tips.
Dimensions: 20.32cm x 6.19cm x 1.9cm
(20Z486)



Capsule Clip-On Holder

Allows quick access to composite during procedure.
Conveniently clips on any side of Tub or B-Lok Tray.
Holds 8 capsules, 4 syringe tubes, and capsule/gun dispenser.
Dimensions: 20.32cm x 5.08cm x 1.42cm
(20Z476)



Trays & Covers



Are you in compliance?

Trays

Our exclusive Safe-Lok® tray cover provides an excellent method of infection control and safety that allows for two benefits other trays do not:

- 1 Prevents any contaminated items from falling off of the used tray during transportation from operator to central sterilization. The contaminated items stay enclosed until you are ready for the cleaning and sterilization process.
- 2 Set up your procedure trays ahead of time with all instrumentation needed. The Safe-Lok® cover securely locks to the tray, not only keeping items clean, but also allowing for safe transportation to chairside. The cover is removed once the patient is in the chair and the procedure begins.

OSHA

Immediately or as soon as possible after use, contaminated reusable sharps shall be placed in puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. (To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.)

CDC

Dirty Instruments must be “contained” prior to transporting them. Containment can be defined as **instrument cassettes** or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.

The Safe-Lok® Tray Cover helps prevent cross-contamination and provide safe transportation back to the sterilization area.



B-Lok Flat Tray

Wide edges for secure grip.
Non-slip feet for stability.
Self-lubricating, high-gloss surface helps prevent staining.
Dimensions: 33.97cm x 24.45cm x 2.22cm

(20Z401) Available in all colors



B-Size Safe-Lok® Tray Cover

Fits B-Lok Flat and Divided Trays only.
Side locking mechanism securely locks cover to tray.
Dimensions: 35.56cm x 25.4cm x 1.91cm

(20Z445)




B-Lok Divided Tray

Holds 12 hand instruments.
Bur holder for 9 FG, RA, or HP burs.
Comes with Insti-Grip.
Non-slip feet for stability.
Self-lubricating, high-gloss surface helps prevent staining.
Holder for Syringe Tip or Saliva Ejector.
Dimensions: 33.97cm x 24.45cm x 2.22cm

(20Z451) Available in all colors



FOR B-LOK DIVIDED TRAY



Insti-Grip

To be used with the B-Lok Divided Tray only.
Holds instruments quietly and safely in place.
Dimensions: 16.83cm x .64cm x .95cm

(20Z454I)



Mini Tray

Self-lubricating, high-gloss surface helps prevent staining.
Non-slip feet for stability.
Wide edges for secure grip.
Mini Trays are great for hygiene/exam setup.
Dimensions: 23.81cm x 16.19cm x 2.22cm

(20Z101) Available in all colors



Mini-Lok Tray Cover

Side locking mechanism securely locks cover to tray.
Dimensions: 25.08cm x 16.83cm x 2.54cm

(20Z103)



Our E-Z Jett® Cassettes fit on Zirc Trays with locking cover! See page 13 for E-Z Jett® Cassettes.



The most efficient method of infection control.

B-Size Tray Cover (Non-Locking)

Dimensions: 35.56cm x 25.72cm x .79cm
(20Z441)



Mini Tray Cover (Non-Locking)

Dimensions: 23.02cm x 17.15cm x 2.22cm
(20Z103)



Multi-Mod Racks

SAVE ON SPACE & Organize Your Office

Tub, tray, and cassette racks are very effective additions to a dental office. They save on counter space and several tubs or trays with covers can fit into one rack. These racks have non-slip feet so they are very stable on any counter or cabinet surface. With a chrome-plated finish, these racks can fit any dental office design.

Multi-Mod Rack 8-Place

Holds 8 Trays or 4 Tubs with covers.
Adjustable to hold Trays, Tubs, and Cassettes.
Sturdy metal construction; chrome-plated.
Non-slip feet.
(21Z410)

S

Multi-Mod Rack 6-Place


Holds 6 Trays or 3 Tubs with covers.
Adjustable to hold Trays, Tubs, and Cassettes.
Sturdy metal construction; chrome-plated.
Non-slip feet.
(21Z105)

S

Multi-Mod Rack (8-Place)



Multi-Mod Rack (6-Place)



Adjustable to fit your needs!



Can adjust to hold Mini Trays and Cassettes



WHY USE Zirc's Cassettes?

A low-cost investment with a quick return! Individual instrument handling takes time away from treating patients and increases overall processing time. Zirc's cassettes will prolong your instrument's life and save you time without breaking the bank.

In a 2014 Implementation of Cassettes Case Report, Marie T. Fluent, DDS, writes how the dental office benefited from the transition from no cassettes to using them with every procedure: "Drs. Dale and Thomas were both grateful that they chose to implement ... an Instrument Management System for the organization, processing, and handling of their instruments within their practice. They agreed that cassettes improved all aspects of instrument management... Initially, their opinions were varied; however, at the conclusion of the one-year anniversary, no clinical staff could imagine processing instruments in any other fashion. Cassette usage had become organized, efficient, safe, and routine." Read the article in its entirety at: www.dentaleconomics.com.

- Designed for superior water flow and cleaning - open top, round holes, and radius corners
- Keeps instrument setups together through sterilization and procedure processes
- Protects instruments from extensive damage and prolongs their life
- Non-slip feet for stability and safety
- AUTOCLAVABLE
- Available in 16 colors for color coding by procedure

Sterilization Pouch Size

Cassette	Pouch Size
E-Z Jett® 5-Place	13.34cm x 30.5cm
E-Z Jett® 8-Place	18.1cm x 35.1cm
E-Z Jett® 8-Place Slim	18.1cm x 35.1cm
E-Z Jett® 10-Place	25.4cm x 38.1cm




Learn more at
www.zirc.com/instrumentcassettes

*Provide pleasing presentation chairside.
Minimize stress for patients.*




**Keep the
workplace
organized and
confusion to a
minimum!**

AUTOCLAVABLE Instrument Holders




Syringe Tip Holder

Syringe Tip Holder 5pk
Dimensions: 1" x 3/8" x 3/8"
(50Z926)




E-Z Jett® Cassette 5-place

Holds 5 instruments.
Great for a basic exam set up!
Syringe Tip Holder included.
Covers stay attached to the cassette and fold back easily.
Fits Statim Sterilizer.
Dimensions: 20.32cm x 9.74cm x 3.02cm
(50Z925) **Available in all colors**



NEW PRODUCT

Size comparison with the 8-Place E-Z Jett® Cassette. See page 3 for more details.




E-Z Jett® Cassette 8-place

Holds 8 instruments.
Covers stay attached to the cassette and fold back easily.
Syringe Tip Holder included.
Dimensions: 20.15cm x 11.26cm x 2.85cm
(50Z927) **Available in all colors**




E-Z Jett® Compact 8-place Slim

Holds 8 instruments.
Covers stay attached to the cassette and fold back easily.
Compact size. Fits the Statim.
Stackable in Ultrasonic.
Syringe Tip Holder included.
Dimensions: 18.09cm x 9.84cm x 1.58cm
(50Z911) **Available in all colors**



Instrument Mat

Holds 10 hand instruments.
Reverse side holds 4 surgical instruments.
Non-slip, soft material.
Small size allows easy access to instruments.
Organizes instruments in order of procedure use.
Can be cut to fit Mini-Tray.
Dimensions: 17.5cm x 5.08cm x 0.95cm
(20Z464) **Available in all colors**




Instrument Strap

Mini Mat

Holds 8 instruments.
Organizes instruments at chairside and with sterilization.
Non-slip feet for stability.
Dimensions: 9.84cm x 8.89cm x 1.27cm
(20Z480) **Available in all colors**

Instrument Strap 2pk

Locking strap secures instruments.
Dimensions: 9.20cm x 1.27cm x .79cm
(50Z917) **Gray Only**



E-Z Jett® Cassette 10-place

Holds 10 instruments.
Accessory compartment in the middle for extra storage.
Covers stay attached to the cassette and fold back easily.
Dimensions: 20.16cm x 17.78cm x 2.85cm
(50Z930) **Available in all colors**



Steri-Containers

AVAILABLE SPRING 2017


Performance Tool Organizer

A deeper, roomier cassette made for large and specialty instruments. Design your own layout with our customizable, removable feet (also sold separately). Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Ideal for:

- Surgical elevators and forceps
- Orthodontic pliers
- Hemostats
- Scissors
- Any instrument that does not fit in your regular cassette



Closed



Open

Performance Tool Kit

Covers stay attached to the cassette and fold back easily.
Dimensions: 20.8cm x 10.8cm x 3.1cm
(50Z940) Available in all colors



Performance Tool Kit Replacement Feet

Two Strips of 10
(50Z944) Gray only




Holds your large and specialty tools

E-Z Jett® Cassettes fit on Zirc Trays with locking cover! Check out our Trays & Safe-Lok® Cover on page 10.



Steri-Container

Holds 10-14 instruments.
Keep instrument setups together through sterilization processes.
Dimensions: 20.63cm x 4.76cm x 4.76cm
(50Z900) Available in all colors



Compact Steri-Container

Holds 8-12 instruments.
Easy open latch, just push down.
Statim 2000 holds 7 Compact Steri-Containers.
Dimensions: 18.09cm x 3.81cm x 3.81cm
(50Z905) Available in all colors



Bur Guards



PROTECT YOUR BURS Steri-Bur Guards




Steri-Bur Guards and Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it also can be positioned at three different angles to give you easy access at chairside.

Steri-Bur Guard

- Fits FG, RA, and HP Burs
- Adjustable height for multiple-length burs (Surgical Bur Guard not adjustable)
- Holds bur in position for maximum cleaning effectiveness
- Drain holes help prevent rusting of burs



12-Hole Bur Guard

Accommodates burs from 12mm-30mm.
Dimensions: 7.3cm x 1cm x 3.5cm
(50Z406) Available in all colors



22-Hole Bur Guard

Accommodates burs from 12mm-30mm.
Dimensions: 12.2cm x 1cm x 3.7cm
(50Z410) Available in all colors



12-Hole Tall Bur Guard

Accommodates burs from 30.9mm-37.1mm.
Dimensions: 7.3cm x .95cm x 4.3cm
Ideal for nose-cone lab burs
(50Z408) Available in all colors

A Short Bur Adapter is included to stop any shorter burs from falling out.

Universal Bur Adapter 5pk

Can be used with 12- or 22-Hole Bur Guard. Enables long and short burs to be used in the same bur guard without short burs falling out. (50Z412)



Endo Guards & Holders



Endo Guards

A unique delivery system to organize, store, and sterilize endodontic files and reamers. The Endo Guard swivels 360° with stops at every 45° for operator's choice. Large slotting provides complete visibility and draining during cleaning.

Steri-Endo Guard

- Holds most lengths for staging, cleaning, and sterilizing files and reamers
- Millimeter gauge measures and positions rubber stopper
- Counter on locking bar tracks the number of sterilizations
- Right- or left-hand use; optional angling
- Available in Hand or Engine styles
- Single-use foam insert cleans/holds files (packaged with 12)



Steri-Endo Guard Engine

Organizes up to 11 engine and 5 hand files/reamers.
Dimensions: 13.65cm x .95cm x 5.56cm
(50Z452) Available in all colors

Steri-Endo Guard Hand

Organizes up to 16 hand files/reamers.
Dimensions: 13.65cm x .95cm x 5.56cm
(50Z450) Available in all colors



Foam Inserts

Single use.
(50Z451) 48pk
(50Z451V) 144pk



Endo Assist

Built-in millimeter gauge measures and positions rubber stopper. For right- or left-hand use. Small and large rings are included to accommodate most finger sizes. Flat bottom enables the Endo Assist to stand freely on its own. Single-use foam insert cleans/holds files (packaged with 12).
Dimensions: 3.81cm x 4.45cm x 4.13cm
(50Z460) Available in all colors

Unique organizer rests on your finger.



Assist Stand

Dense foam cleans bio-debris faster and more easily. Holds Endo Instruments in an upright position. Square shape fits great in procedure tub. Thicker foam prevents solution from leaking. Rounded corners for easy cleaning.
Dimensions: 4.76cm x 4.76cm x 3.16cm
(50Z462) Black Only



Assist Stand Foam

Single use.
(50Z466) 48pk
Dimensions: 26.67cm x 19.685cm x 2.22cm
(50Z466V) 144pk
Dimensions: 26.67cm x 19.685cm x 2.22cm



Bur Blocks



Magnetic Bur Blocks

- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- Non-slip feet for stability
- Drain holes help prevent rusting of burs



8-Hole Bur Block

Covers sold separately.
 Dimensions: 3.81cm x 3.81cm x 1.27cm
 (50Z400) Available in all colors



14-Hole Bur Block

Covers sold separately.
 Dimensions: 7.46cm x 3.81cm x 1.27cm
 (50Z403) Available in all colors



Bur Block Covers

- Covers available to maintain cleanliness of burs during storage
- 1" covers for FG and RA burs; 2" cover (14-hole only) for lab and surgical burs
- Covers are not autoclavable



8-Hole Bur Block Covers 1" High - Standard

Dimensions: 4.13cm x 4.13cm x 2.54cm
 (50Z401)



14-Hole Bur Block Covers 1" High - Standard


Dimensions: 7.94cm x 4.13cm x 2.54cm
 (50Z404)



14-Hole Bur Block Covers 2" High - For Surgical Burs

Dimensions: 7.94cm x 4.13cm x 5.72cm
 (50Z405)






How to effectively identify your instrument setup.

- 1 color for practitioner:**
Neon Orange shows ownership to the practitioner.
- 1 color for procedure:**
Neon Blue is for identifying the procedure.

Did you know?
The diagonal setup of the orange Rings shows the order of how the instruments should be used. This eliminates confusion between the doctor and assistant.

Place your instrument to determine which ring size you need.



Now available in **XL size!**



E-Z ID Rings

- Prepackaged in an easy-to-use dispenser
- Available in Small, Large, and XL sizes to fit all instrument handles

25pk

Available in all colors.

- (70Z100)
- (70Z200)
- (70Z000)

Jewel System 200pk

Assorted colors (E, J, T, I, M, and U).

- (70Z115)
- (70Z215)
- (70Z015)

Vibrant System 200pk

Assorted colors (N, O, P, Q, R, and S).

- (70Z110)
- (70Z210)
- (70Z010)

Classic System 200pk

Assorted colors (A, B, D, G, and I).

- (70Z105)
- (70Z205)
- (70Z005)



E-Z ID Tape

- Tape dispenser has built-in tape cutter.
- Dispenser roller prevents tape from twisting.



E-Z ID Tape Roll

Available in all colors.
3m roll.
(70Z300)



E-Z ID Tape System - Jewel

Assorted colors (E, J, T, I, M, and U).
Includes : 0.91m rolls.
(70Z315)

E-Z ID Tape System - Vibrant

Assorted colors (N, O, P, Q, R, and S).
Includes : 0.91m rolls.
(70Z310)

E-Z ID Tape System - Classic

Assorted colors (A, B, D, G, and I).
Includes : 0.91m rolls.
(70Z305)

Label Containers!
Keep containers organized
by labeling them per content.



Available in U.S.A. only

ZOBE

the experience


Zirc's **O**rganizational **B**ox of **E**fficiency

Why Z.O.B.E.?

- The ZOBE kit ships to your practice at NO cost – Free ZOBE, free shipping, free return
- ZOBE contains over 25 of our top-selling products for organization
- Have fun exploring the products – **Please Note: Z.O.B.E. is for demonstration purposes ONLY**
- Your ZIRC efficiency expert consults with you and listens to your needs
- Your ZIRC expert reviews the amount of patient flow for each procedure
- You receive a custom plan that's right for you
- Virtual consultations available

Visit us at www.zirc.com/zobe





Dreaming of a stress-free, fun, and organized practice?
Let us send you Z.O.B.E.

The ZOBE kit ships to your practice at NO cost.* Free ZOBE. Free Shipping. Free Return.

*Available in USA only, contact ZIRC for more information.



Pedo/Small Adult



Medium/Large



The Efficiency of Hands-Free Isolation

A hands-free isolation device that retracts, isolates, and evacuates! Work continuously without the repeated use of suction, minimize your isolation inventory, and gain a quick return. With two convenient color-coded sizes, these one-time use isolation devices are the perfect cost-effective tool for your dental office.

- Built-in bite block for patient comfort
- Connects directly to HVE Valve
- Keeps work environment dry

Pre-Assembled

The new mr. thirsty® one-step is ready to use and is single use! Just insert into your HVE Valve (one-step!).

2 Sizes Available

mr. thirsty® one-step comes in two convenient sizes, Pedo/Small Adult and Medium/Large. It can also be trimmed for a custom fit!



10pk Trial Kit

10mm Adaptor (50Z98710)
16mm Adaptor (50Z98716)
Includes 5 of each size:

- (5) Pedo/Small Adult
- (5) Medium/Large



25pk

- Pedo/Small Adult
10mm Adaptor (50Z986SM10)
16mm Adaptor (50Z986SM16)
- Medium/Large
10mm Adaptor (50Z986LG10)
16mm Adaptor (50Z986LG16)

ALSO AVAILABLE IN 100 PACK



Comfort Kit

The Comfort Kit's lightweight tubing allows your heavy HVE Valve/Tubing to stay mounted to the equipment.

Includes:
Rotary plastic HVE Valve
1.83m lightweight tubing
mr. thirsty® one-step NOT included
(50Z992)



Take comfort to the next level with lightweight tubing!



mr. thirsty® one-step benefits everyone! Here's how...



DOCTOR

Doctor
25% average decrease in procedure time. Hooks directly to HVE Valve. Single-Use.



HYGIENIST

Hygienist
Reduce chair time by as much as 20%. Hands-free continuous suction. Ease patient anxiety.



ASSISTANT

Assistant
Free up time. See another patient in the next room. Catch up with sterilization.



PATIENT

Patient
Improved patient comfort. Non-threatening appearance and improved patient safety.




REDUCE CHAIR TIME BY AS MUCH AS 20%




How does mr. thirsty® one-step stack up against the competitors?

mr. thirsty® one-step Competitor Brand

Cost-Effective	✓	
Large Upfront Investment		✓
No Maintenance or Sterilization Required	✓	
Can be Trimmed to Fit	✓	
Colored Product for Oral Cavity Contrast	✓	






To read a full product evaluation visit...
<http://www.dentalproductshopper.com/relaxed-fit-insti-dam>

**Easy to use, non-threatening,
 and comfortable for patient.**

Insti-Dam®

Features for BOTH Relaxed Fit and Original Insti-Dam®

- Built-in flexible medium gauge material that comes pre-assembled with frame and pre-punched hole
- Pre-punched hole makes positioning easy and allows eyes and nose to be open and comfortable
- Helps eliminate tearing (additional holes can be punched)
- Compact design fits outside patients lips, easy to use, non-threatening and comfortable for patient
- Radiographs may be taken by folding it to the side
- Single use



Latex-Free Insti-Dam® Relaxed Fit 20pk

Ideal for posterior isolation.
 The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing, or pulling.
 Dimensions: 10.79cm diameter

Patent #7044738
 (50Z457)

Blue Latex-Free Insti-Dam® 20pk

Ideal for anterior isolation.
 Made with polyisoprene.
 Dimensions: 10.79cm diameter
 (50Z459)

Natural Latex Insti-Dam® 20pk

Ideal for anterior isolation.
 Dimensions: 10.79cm diameter
 (50Z455)

Lingua-Fix® & SE Cushion



Lingua-Fix® Saliva Ejector

Lingua-Fix® provides a dry work area.
Holds tongue in a steady position to protect from injuries.
Simple operation... easy to place and remove.
Comfortable, no sharp edges or corners.
Large suction surface avoids blockage and is more comfortable.
Reduces sedimentation of drilling by-products in pipes and tubing.
Single use only.
Comes pre-assembled.
(50Z953S) **10pk**
(50Z953) **50pk**



SE Cushion

Used in conjunction with the Saliva Ejector.
Prevents tissue grabbing and tissue irritation.
Increases patient comfort.
Improves range of suction.
Quiets the noise of the suction.
Single use.
(50Z970) **100pk**



Can be
used with
all Dental
Dams!

**For comfortable and efficient
suction or evacuation.**

The Relax Factor:

The extra "relaxed" material in the Relaxed Fit Insti-Dam® makes for quick posterior placement - no stretching, tearing or pulling.



Insti-Dam® Dispenser

Made for Insti-Dam® Original ONLY
Can be mounted to wall or cabinet or can sit on a counter.
Available in White (A) for latex or Neon Blue (N) for latex-free.
Non-slip rubber bottom for stability.
Holds 35 Insti-Dams®.
Dimensions: 12.38cm x 12.38cm x 15.56cm
(50Z471) **Color Code: A & N**

AM S



A TV without HD is like a...
A MIRROR WITHOUT HD



Crystal HD


Experience high-definition quality. The brighter the workspace, the better the illumination and visibility. With ZIRC's Crystal HD® lens, you are visually working in a 40% brighter environment compared to the Rhodium lens that is currently the industry standard. If you haven't made the change to HD - what are you waiting for?

**40%
BRIGHTER**
Than Rhodium
Mirrors



Your mirror **VS** Crystal HD®

CUTTING EDGE IMAGERY



Rhodium = Reflects color 75% accurately
Crystal HD® = Reflects color 99% accurately

Provides exceptionally bright,
shadow-free images!



Crystal HD®

Available in multiple styles!

Designed to reduce eye strain and fatigue.

Autoclavable



As a microscope user, I can understand the importance of clear vision. ZIRC mirrors significantly improved the brightness of my image, giving my eyes highest quality vision. It's true - brightness makes all the difference.

- Dr. Dimitris Tsanaksidis, Microscopic & Aesthetic Dentistry

CRYSTAL HD® MIRROR + MIRROR MAGIC®
The Perfect Pair

Pair our Crystal HD® lens with our mirror magic® System for the highest quality vision and best procedure outcome. See page 30.




Crystal HD® Mouth Mirrors & Handles

Flush surface, no debris build-up can occur



Resin

Gap where debris build-up occurs



Stainless Steel

	RESIN	STAINLESS
Autoclavable	X	X
Antimicrobial protection	X	
Lens flush with mirror frame	X	
No galvanic shock on patient	X	
Rust-free	X	

Why use Zirc's Resin Mirrors?



Assorted Soft Grip 12pk

Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

#4 (50Z363-NEON)

#5 (50Z364-NEON)

Jewel Assorted

2 of each color: E, J, T, I, M, and U.

#4 (50Z363-JWL)

#5 (50Z364-JWL)



Assorted Thin Grip 12pk

Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

#4 (50Z353-NEON)

#5 (50Z354-NEON)

Jewel Assorted

2 of each color: E, J, T, I, M, and U.

#4 (50Z353-JWL)

#5 (50Z354-JWL)



Black/Gray Soft Grip 12pk

#4 (50Z363)

#5 (50Z364)

Black/Gray Soft Grip 4pk

#4 (50Z361)

#5 (50Z362)



Handles



Ergo Grip Handle

(50Z347) Color Code: B-S, and K (black)





Resin Screw-on Heads 12pk

Flush surface, no debris build-up can occur.
Simple stem.

1-Color 12pk

Available in all colors and K (black).

#4 (50Z392)

#5 (50Z393)

Vibrant Assorted

2 of each color: N, O, P, Q, R, and S.

#4 (50Z392-NEON)

#5 (50Z393-NEON)

Jewel Assorted

2 of each color: E, J, T, I, M, and U.

#4 (50Z392-JWL)

#5 (50Z393-JWL)

Classic Assorted

A, B, D, G, I, and K (black).

#4 (50Z392-PSTL)

#5 (50Z393-PSTL)



Vibrant Assorted Double-Sided Mirror Head 4pk

1 of each color: N, O, R, and S.
Simple stem.

#4 (50Z394-NEON)



Gray Double-Sided Mirror Head 4pk

Simple stem.

#4 (50Z394)



Stainless Steel Head 12pk

Simple stem.

#3 (50R389)

#4 (50R390)

#5 (50R391)



Stainless Steel Head 1pk

Simple stem.

#0 (50R388)



Zirc's Mirrors & Handles are autoclavable!

mirror magic® Anti-Fog Solution



without mirror magic®



with mirror magic®



mirror magic®
anti-fog system

Anti-Fog Solution removes debris from dental instruments!



mirror magic® Trial Kit

2-6ml bottles of Anti-Fog Solution
24 Swipe Pads
(5OR200X)

Increase Your Chairside Efficiency

The mirror magic® Anti-Fog System prevents fogging on mouth mirrors and offers the practitioner a quick method for removing debris during procedures. The disposable Swipe Pads stick to the back of a gloved hand or a patient bib for easy access. The Anti-Fog System requires only a small amount of solution and is economic in size, making it easy to store.



Watch the magic in action
at www.zirc.com/mirrormagic

Combine mirror magic® with Crystal HD® Mouth Mirror for the best view in dentistry.





"The mirror magic® Anti-Fog System decreases fogging while also decreasing the build-up of debris on my mirror. This affords me outstanding visualization, which is crucial to performing quality dentistry. mirror magic® has truly made a difference in my practice and is one of those features I didn't realize I needed; and now I can't imagine practicing without it."

— Andrew M. Goldsmith, DDS, DICOI, FIALD



10pk Anti-Fog Solution (6ml ea.)

1 bottle lasts approx. 50 patients.
(50R205)



Swipe Pads 96pk

Single use.
(50R210)



mirror magic® Organizer

Holds: 48 Swipe Pads and 1 - 6ml bottle of Anti-Fog Solution.
Store on operatory counter or in drawer.
Non-slip bottom for secure placement.
Note: box only, does not include mirror magic.
Dimensions: 21.27cm x 5.87cm x 2.86cm
(50R201)



Clarity Collection

Includes:
Anti-Fog Solution 10pk
Swipe Pads 96pk
1 Crystal HD® Mouth Mirror
1 mirror magic® Organizer
(50R220)

The Perfect Pair

The superior visual acuity and 40% brighter Crystal HD® lens combined with the quick debris removal and Anti-Fog Solution of mirror magic® have proven to be the perfect pair. Increase chairside efficiency and speed up your process with this dynamic duo.



Denture/Retainer Boxes



Retainer Box

One piece construction.
Solid locking mechanism.

Retainer Box 2.54cm Deep 12 pk

Dimensions: 7.62cm long x 2.54cm deep

(25R500) 12pk

Color Code: J, N, O, P, R, S, Assorted

Retainer Box 3.81cm Deep 12 pk

Dimensions: 7.62cm long x 3.81cm deep

(25R550) 12pk

Color Code: J, N, O, P, R, S, Assorted

S

Solid locking mechanism so your retainer stays put!

Great colors to choose from!

Denture Box

One piece construction.
Solid locking mechanism.

Denture Box 12 pk

Dimensions: 9.53cm long x 4.45cm deep


(30R800) 12pk

Color Code: B, C, G, H, J, R, Assorted

S



Saliva Ejector Accessories & MEGA Products



Rotary

Lever

Hands free isolation!
See page 25.

Saliva Ejector Valve (with Shut-Off)



Quick disconnect with swivel.
Rubber tip with protective screen.
Provides variable control from full on to full off.

Saliva Ejector Valve with Rotary On/Off Control

Dimensions: 7.62cm x 1.59cm x 1.91cm
(23E360) **Gray Only**

Saliva Ejector Valve with Lever On/Off Control

Dimensions: 7.62cm x 1.59cm x 1.75cm
(23E365) **Gray Only**

Saliva Ejector Valve Lever Only

Dimensions: 1.59cm x .48cm x 1.59cm
(23E364) **Gray Only**

Saliva Ejector Replacement Tip

Dimensions: 1.59cm x 1.59cm x 2.54cm
(23E363) **Gray Only**



Saliva Ejector Screens

Single use.

Saliva Ejector Screens

Dimensions: .95cm Diam. x .95cm
(50Z952S) **25pk**
(50Z952P) **250pk**



MEGA Products

MEGA Gracey

39.37cm in length.
Great teaching tool.
A gift for any dental personnel.
(50Z312) **Mustard Only**



MEGA Mirror

(Patient hand mirror).
6.35cm lens and 39.37cm in length.
Keeps children entertained.
Great gift for any dental personnel.
(50Z309) **Color Code: N, P, Q, R, S**



Fun gift ideas for any dental personnel.
Keeps children entertained.




Water Wise & Steri-Soaker



Water Wise Starter Kit

Holds 1 liter.
Includes bottle, adaptor, lid, neon yellow and neon blue bubbles (water level indicator), and tube.
(23Z100)



Water Wise Replacement Bottle

Holds 1 liter.
Includes bottle, adaptor, lid, neon yellow and neon blue bubbles (water level indicator), and tube.
(23Z105)



Water Wise vs. Traditional Water Bottle

	Water Wise	Water Bottle
Easy to clean	X	
Antimicrobial	X	
Easy-to-fill large mouth	X	
Quick release	X	
Level indicator color coded for: Sterile and distilled water, Medicament	X	
Strong	X	X
Easy to grip	X	
Sealable lid	X	
No threads to strip	X	
Pressure relief valve	X	

We strive to make things easier around the dental office.

The Steri-Soaker

An efficient, high-volume germicide tray designed for maximum use of chemicals to clean/disinfect instruments. It is designed for easy cleaning and is large enough to accommodate ZIRC Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.

Steri-Soaker

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters. Hinged fitted cover.

Safe-to-use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines.

Dimensions: 35.24cm x 24.45cm x 13.02cm

(50Z805) Color Code: B, I



Designed for easy cleaning!



E-Z Access Shelf/Accessories


E-Z Access Shelf

E-Z mounting to standard equipment post (4.45cm-6.98cm); No tools needed.
 Will not damage post's finish.
 E-Z adjustment to multiple positions.
 Contoured edge to contain spills and materials.
 Stain resistant finish; surface disinfect only.
 Weight capacity of 11.33kg
 Workable tray area: 34.93cm L x 26.04cm W.
 Standard E-Z Access Shelf as is accommodates 4.45cm-5.72cm post.
 Rotates 360°.

(20Z420) Color Code: A, G, I



Provides easy access work space for auxiliary items.



E-Z Access Shelf Adapters

E-Z Access Shelf Adapter for Smaller Posts

Adapts to posts 3.49cm-6.99cm
 Dimensions: 5.13cm x 4.76cm x 1.91cm
 (20Z422) Color Code: A, G, I

3-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 5.72cm-6.03cm
 Dimensions: 10.16cm x 5.08cm x 5.08cm
 (20Z437) Color Code: A, G, I

4-1/2" Bolt Adapter Kit for Larger Posts

Accommodates post sizes from 6.03cm-6.99cm
 Dimensions: 11.43cm x 5.08cm x 5.08cm
 (20Z438) Color Code: A, G, I



Cup Holder-Snap On

Fits all standard cups.
 No tools needed.
 Dimensions: 7.94cm x 6.35cm x 2.54cm
 (20Z436I)




E-Z Access Shelf Surface Mount

Mount on horizontal flat surface only.
 Dimensions: 9.53cm x 9.53cm x 7.30cm
 (20Z424) Color Code: A, G, I



Our products are able to fit every need!



Wall or Cabinet Mount

Allows you to mount the E-Z Access Shelf onto a wall or the side of a cabinet. 90° angle of mount permits 180° movement of shelf on cabinet or next to wall.
 Dimensions: 16.51cm x 14.61cm x 12.38cm
 (20Z426) Color Code: A, G, I



8" Extension Arm

Enables shelf to extend over patient.
 Shelf has strength to hold 25 lbs. when connected to arm.
 May be used with Wall Mount, Surface Mount, or directly to post.
 Dimensions: 25.4cm x 10.16cm x 11.43cm
 (20Z428) Color Code: A, G, I



Cabinet Trays & Lab Pans

AM S 

Cabinet Trays



No 16A – Hand Instruments (Regular Well)

Dimensions: 20cm x 9.53cm x 2.38cm
(20Z201) Color Code: A, B




No 16 – Hand Instruments (Deep Well)

Dimensions: 19.69cm x 9.53cm x 2.54cm
(20Z202) Color Code: A, B



No 17 – Clamps, Retainers, etc.

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z203) Color Code: A, B



No 18 – Finishing Strips

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z204) Color Code: A, B



No 19 – Pliers, Scissors, etc.

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z205) Color Code: A, B



No 20 – Discs, Wheels, Stones

Dimensions: 20cm x 9.53cm x 1.91cm
(20Z206) Color Code: A, B



Wall-Hanging

AM S 

Lab Pans

Organizes and transports impressions, models, dentures, and crown and bridge appliances.
Stainless-steel spring clip holder.
Color-coded.

Lab Pan

Dimensions: 19.05cm x 14.61cm x 17.78cm
(17Z101) Color Code: A, B, D, G, M

Wall-Hanging Lab Pan

Same as standard version but with strip on top back for hanging on any vertical surface.

Saves on installing expensive shelving.
Dimensions: 19.69cm x 14.61cm x 17.78cm

(17Z102) Color Code: A, B, D, G, M



Wall Strip Lab Pan

Drilled and countersunk with screws and wall anchors.
45.72cm Wall Strip.


Holds 3 Lab Pans.

Dimensions: 45.72cm

(17Z104)

S

Save on installing expensive shelving.



Deluxe Cotton Roll Dispenser

Heavy-duty no-skid base.
Dispenses up to three rolls at a time.
Extra sturdy construction.
Dimensions: 10.64cm x 5.24cm x 10.16cm
(28R800) **Color Code: A, B**

S



Cotton Roll Holder

Stainless pins embedded in weighted bottom.
Neatly holds cotton rolls in an upright position.
Hinged lid helps keep cotton rolls clean, prevents spilling.
Dimensions: 9.21cm x 8.26cm x 6.03cm
(28R900) **Color Code: A, B**

S



Pellet Dispenser

Spring loaded – pushes pellets to the surface.

Dispenser for Small Pellets (sizes 3, 4, 5)
Dimensions: 4.76cm x 4.76cm x 5.24cm
(28R825) **White Only**

Dispenser for Medium Pellets (sizes 0, 1, 2)
Dimensions: 4.76cm x 4.76cm x 5.24cm
(28R850) **White Only**

S

Sponge Dispensers

Spring action releases pads one at a time.
Bottoms are weighted to prevent tipping.

5.08cm x 5.08cm Sponge Dispenser

Dimensions: 6.35cm x 6.35cm x 7.3cm
(28R950) **Color Code: A, B**

10.16cm x 10.16cm Sponge Dispenser

Dimensions: 11.43cm x 11.43cm x 7.3cm
(28R975) **Color Code: A, B**

S



Matrix Dispenser

One roll of 1cm x 11m Matrix clear tape 0.002 gauge.
Allows Matrix to be contoured (u-shaped) as it dispenses.
Stainless steel cutter allows for desired length.
Disposable.
Dimensions: 5.08cm x 5.39cm
(55Z200) **Color Code: A, B**

S

E-Z Storage/Handi-Hopper

E-Z Storage Tub System

Hinges allow cover to stay open on its own!
Each Tub comes with one vertical divider and one half-depth adjustment divider.
Hinged cover comes in white or clear.

E-Z Storage Organizer (with white cover)

Dimensions: 14.92cm x 12.38cm x 14.61cm
Rack sold separately.
(29R605)

E-Z Storage Organizer (with clear cover)

Dimensions: 14.92cm x 12.38cm x 14.61cm
Rack sold separately.
(29R610)

Replacement Hinges (pair)

Dimensions: 2.22cm x 1.27cm x .32cm
(29R627)

S



E-Z Storage Single Rack

Holds storage organizer Tub.
Tub can be placed flat or at an angle for easy access of supplies.
Wall or cabinet mounted with optional mounting bracket.
Sturdy metal construction; plastic-coated.
White, non-slip feet.
Dimensions: 17.78cm x 12.7cm x 12.06cm
(29R650)

S

E-Z Storage Mounting Bracket

White metal bracket for the E-Z Storage Tub single rack.
Pre-drilled holes for cabinet or wall mount.
Dimensions: 17.78cm x 12.7cm x 12.07cm
(29R651)

S



AM S

Mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste.

Handi-Hopper (with Mounting Bracket)

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean.
Slotted mounting bracket (included) for easy removal and replacement.
To use for infectious waste, order red (Biohazard) liners.
To use for other waste, order clear liners.

Handi-Hopper

Dimensions: 17.78cm x 10.16cm x 16.19cm
(20Z404) Color Code: G, I

Handi-Hopper Bracket

Dimensions: 5.72cm x .64cm x 2.54cm
(20Z408) Color Code: G, I



Handi-Hopper Liners (Single-Use)

Handi-Hopper Liners 100/Box (Red Biohazard)
(20Z405) 100pk

Handi-Hopper Liners 100/Box (Clear)
(20Z410) 100pk

AM S



Waste-Drop Ring

Waste-Drop Ring for cabinet top.
Fits 12.7cm hole.
Made of self-lubricating, easy-to-clean resin.
Dimensions: 14.29cm x 14.29cm x 7.62cm
(22Z100) Color Code: A, G, I

Cuspidor/Dispos-A-Screen Traps & Needle Capper



Reduce possibilities of infection and contamination.

Cuspidor Screens AM

Single use.

Cuspidor Screens Large Universal – 6.35cm (gray)

Adec, Ampco, Belmont, Dab Atlante, Dansereau, DCI Den-Tal-Ez, Forest, Knight, Marus, MDT, Proma, Ritter, Schein.

(51Z205S) **Gray 25pk**

(51Z205V) **Gray 144pk**

Cuspidor Screens Small – 4.13cm (white)

Adec, Cascade and Performer cuspidor bowls.

(51Z200S) **White 25pk**

(51Z200V) **White 144pk**

Dispos-A-Screen Traps S

Avoid contact with potentially hazardous material. Improves performance of suction system. Color coded by size. Single use.

(25pk)

(144pk)

Blue – 4.5cm

Adec models with black plastic canisters, plus newer models by Ampco, Belmont, DCI, Forest, Marus MDT, Schein.

25pk (51Z045S), **144pk** (51Z045V)

Green – 4.76cm

Ampco, Belmont, Biotec, Chayes, Dabi-Alante, Dansereau, DCI, Den-Tal-EZ, Dentech, Forest, Knight, Pelton & Crane, Proma and others.


25pk (51Z123S), **144pk** (51Z123V)

Gray – 7cm

Midwest 210 Models.

25pk (51Z070S), **144pk** (51Z070V)

Left to right: Blue, Green, Gray, Yellow, Beige, Mauve and Coral



Yellow – 5.72cm

Adec, model with metal canisters.

25pk (51Z040S), **144pk** (51Z040V)

Beige – 7.3cm

Den-Tal-EZ models with metal canister.

25pk (51Z060S), **144pk** (51Z060V)

Mauve – 3.18cm

Adec, Cascade, and Performer.

25pk (51Z080S), **144pk** (51Z080V)

Coral – 6.67cm

Pelton and Crane Spirit Unit with Plastic Canister.

25pk (51Z050S), **144pk** (51Z050V)



Needle Capper

One-handed capping and recapping.

Syringe is for use; at convenient position in Needle Capper.

Base of device is secured to tabletop or tray with adhesive strip (included).

Needle is uncovered at same time syringe is withdrawn from Capper.

After use, syringe is replaced for recapping.

Needle Capper available in Standard and Long sizes.

Either size tube will fit standard base.

Standard Needle Capper

(base, one tube, adhesive tape)

Dimensions: 5.56cm x 3.18cm x 3.18cm

(27R700) **White Only**

Long Needle Capper

(base, one tube, adhesive tape)

Dimensions: 6.35cm x 3.18cm x 4.45cm

(27R704) **White Only**



Mighty Mixer Bowls, Spatulas & Organizers

Single-use
disposable mixing
bowls.



Mighty Mixer Bowls 36pk (White Only)

Disposable, eliminating cleanup time.
Fits securely into mixing bowl without slippage.
Saves wear on mixing bowl.
Free-standing for use with or without mixing bowl.
Capable for use with electric Alginate mixers.
Dimensions: 12.38cm Diam. x 7.93cm
(50Z505) 36pk



Mighty Mixer Bowls 36pk (Assorted)

Same features as the white Mighty Mixer Bowls.
Dimensions: 12.38cm Diam. x 7.93cm
9 of each color: J, N, P, and R
(50Z505-ASTD) 36pk



Mighty Mixer Spatula 4pk (Assorted)

Same features as the white Mighty Mixer Spatula.
Dimensions: 18.57cm .48cm x 3.18cm
1 of each color: J, N, P, and R
(50Z501-ASTD) 4pk



Mixing Sticks 100pk

Single use.
Dimensions: 11.43cm long
(50Z521) 100pk



Cement Spatula

White plastic, double-end spatula.
For fillings, cement, and pastes.
Self-lubricating plastic for easy cleaning.
Dimensions: 16.19cm x .64cm x .95cm
(50Z502)



Delrin Spatula 10pk

Flexible flat blade.
Non-stick.
Resistant to chemicals.
Great for mixing cements, acrylics, etc.
Dimensions: 11.74cm x .47cm Diam.
(50Z523) 10pk



Mighty Mixer Spatula

Self-lubricating for easy cleaning.
White plastic spatula.
Broad, medium flex blade.
Ideal for alginates and plaster.
Dimensions: 18.57cm .48cm x 3.18cm
(50Z501) White Only



Implant Organizer

Quick and easy recognition. Individual compartments for implants.
Shape of the organizer matches the arch of the mouth for implant placement.
Separate compartments prevent implant appliances from mixing together.
Lubricated plastic enables tooth number to be written in the designated area.
Cover helps keep materials secure.
Dimensions: 17.46cm x 11.75cm x 1.59cm
(20Z490B) Blue Only



No more
confusion!



Restorative Organizer

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays.
Disposable Mixing Wells included for bonding agent and cement.
Accessory area for brushes and materials.
Lubricated plastic enables tooth number to be written in the designated area.
Light protective cover can be used as a light shield and to prevent materials from setting early.
Includes 25pk of Mixing Wells.
Dimensions: 17.46cm x 11.75cm x 1.59cm
(20Z495R) Neon Purple Only



Mixing Well (Single-Use) 144pk
(20Z496V) 144pk



Contact us today to get started!
+1 763.682.6636 | sales@zirc.com